



mocom[®]
Sistemi innovativi di sterilizzazione

Quality System certified by



**RAPID STEAM
STERILIZER**

millennium[®] B_μ

OPERATING MANUAL

UPDATING

The following table lists the updating history of the Operating Manual.
The field "Description" shortly explains the modifications.

Ed.	Rev.	Date	Description
1	0	03-04	First issue

TABLE OF CONTENTS

1. FOREWORDS	1
APPLICABLE EUROPEAN DIRECTIVES	1
USE IDENTIFICATION	1
SCOPE OF THE MANUAL	2
GENERAL WARNINGS	2
2. PACKING CONTENT	3
DIMENSIONS AND WEIGHT	3
CONTENT	3
LIFTING AND TRANSPORT	4
3. GENERAL DESCRIPTION	5
INTRODUCTION	5
GENERAL FEATURES	5
FRONT DEVICES	6
REAR DEVICES	6
COMMAND PANEL	7
LCD DISPLAY	7
OPERATING CYCLE EXAMPLE	8
4. INSTALLATION	9
GENERAL	9
DIMENSIONS FOR BUILD MOUNTING	9
GENERAL INSTALLATION DIRECTIONS	10
MAINS CONNECTION	10
EXTERNAL PRINTER CONNECTION (OPTION)	10

5. FIRST START-UP	11
SWITCHING ON.....	11
INITIAL SELF-TEST	11
ACQUISITION AND AUTOMATIC UPDATING OF THE ENVIRONMENT PRESSURE VALUE	11
STAND-BY STATUS	12
FILLING THE DISTILLED WATER.....	13
RECOVERY TANK MAX WATER LEVEL SIGNALING.....	13
6. SETTING THE EQUIPMENT	14
INTRODUCTION	14
STARTING THE SETUP PROGRAM.....	14
SETUP MODE KEYS FUNCTION.....	14
DESCRIPTION OF THE MENU OPTIONS	16
DEFAULT SETTINGS	18
SETTING THE OPTIONS IN SETUP MODE.....	18
Language	18
Date setting	18
Time setting	19
Setting the preset sterilization Programmes	20
Setting the STAND-BY modes	24
Setting the Printer Options	25
Setting the tank Filling mode	27
Setting the Draining mode	27
Acquisition of the environment pressure value	27
Adjusting the LCD contrast	28
EXIT THE SETUP MODE.....	28
7. PREPARING THE MATERIAL TO BE STERILIZED	29
INTRODUCTION	29
HANDLING THE MATERIAL BEFORE STERILIZING	29
ARRANGEMENT OF THE LOAD	30
8. SELECTING THE STERILIZATION PROGRAM	32
INTRODUCTION	32
ABOUT THE SELECTION.....	32
9. RUNNING THE STERILIZATION PROGRAM	34
GENERAL	34
STARTING THE PROGRAM.....	34
SEQUENCE OF THE PROCESS	35
MANUAL INTERRUPTION OF THE PROGRAM	39
RESULT OF THE PROGRAM.....	40
PRINTING THE DATA REPORT	40
10. PRESERVING THE STERILIZED MATERIAL.....	41
GENERAL	41
HANDLING	41
STORING	41
11. TEST PROGRAMMES	42
OVERVIEW	42

HELIX/BD TEST	42
VACUUM TEST	43
APPENDIX A – TECHNICAL CHARACTERISTICS	46
OVERVIEW TABLE	46
SAFETY DEVICES	47
CHARACTERISTICS OF THE FEEDING WATER	48
APPENDIX B – PROGRAMMES	49
INTRODUCTION	49
OVERVIEW OF THE AVAILABLE PROGRAMMES	50
DIAGRAMS OF THE STERILIZATION PROGRAMMES	51
DIAGRAMS OF THE TEST PROGRAMMES	56
EXAMPLES OF PRINTING REPORTS	57
APPENDIX C - MAINTENANCE	59
GENERAL	59
ORDINARY SCHEDULED MAINTENANCE	59
MAINTENANCE ACTIVITIES	60
Cleaning the gasket and the inner part of the door	60
Cleaning the external surfaces	60
Cleaning the sterilization chamber and accessories	60
Disinfecting the external surfaces	60
Maintenance of the safety valve	60
Cleaning (or replacing) the draining filter	61
Sterilizing the bacteriological filter	61
Replacing the bacteriological filter	61
PERIODIC VALIDATION OF THE STERILIZER	61
APPENDIX D - GENERAL PROBLEMS	62
OVERVIEW	62
ANALYSIS AND RESOLUTION OF THE PROBLEMS	62
APPENDIX E – ALARM SIGNALING	65
OVERVIEW	65
ALARM PROCEDURE	65
Alarm occurring during the sterilization cycle	65
Alarm occurring out the sterilization cycle	66
RESET OF THE SYSTEM	67
LIST OF THE ALARM CODES	68
ANALYSIS AND RESOLUTION OF THE PROBLEMS	70
APPENDIX F – DRAWINGS	77
ELECTRICAL DRAWING	77
HYDRAULIC DRAWING	78
APPENDIX G – DECLARATION OF CONFORMITY	79
APPENDIX H –OPERATOR'S NOTES	80
APPENDIX Z – CUSTOMER SERVICE	81

FOREWORDS

On thanking for the preference granted, M.O.COM. Ltd. Co. hopes that the performances of this product can result of Your complete satisfaction.

In this manual you will find all the procedures for the correct use and the indications for the complete exploitation of the equipment performances.

We remain at your disposal for any more explanation, as well as suggestions turned to improve the product and the service.

Symbols used through the manual



PAY PARTICULAR ATTENTION TO THE PARAGRAPHS MARKED BY THE FINGER SYMBOL.



ATTENTION! THIS SYMBOL POINTS OUT A POTENTIAL DANGER FOR PEOPLE. PLEASE OPERATE ACCORDING TO THE SUITABLE PROCEDURES OF THE MANUAL IN ORDER TO PREVENT POSSIBLE DAMAGES TO THE USER AND/OR TO THIRD PARTIES.



CAUTION! THIS SYMBOL POINTS OUT A POTENTIAL DANGER FOR PROPERTY. ACT ACCORDING TO THE SUITABLE PROCEDURES OF THE MANUAL IN ORDER TO PREVENT POSSIBLE DAMAGES TO MATERIAL, EQUIPMENT AND/OR PROPERTIES.



CAUTION! THIS SYMBOL POINTS OUT A POTENTIAL DANGER DUE TO PRESENCE OF **HIGH TEMPERATURE**.

APPLICABLE EUROPEAN DIRECTIVES

The product described in this manual is manufactured in accordance with the highest safety standards and doesn't represent any danger for the operator if used according to the following instructions.

The product is in accordance with the following European Directive as applicable:

- 73/23/CEE,** for the approximation to the legislation of the Members States related to low voltage equipment (and following modifications).
- 89/336/CEE,** for the approximation to the legislation of the Members States related to the electromagnetic compatibility (and following modifications);
- 93/42/CEE,** concerning the medical devices (and following modifications).

USE IDENTIFICATION

The product described in this manual is exclusively intended for the sterilization of re-usable surgical instruments and material.



THE USE OF THE EQUIPMENT IS STRICTLY LIMITED TO QUALIFIED PERSONEL. NO GROUNDS JUSTIFY THE EQUIPMENT USE OR HANDLING BY UNSKILLED AND/OR UNAUTHORISED PERSONEL.

THE DEVICE MUST NOT BE USED FOR THE STERILIZATION OF FLUIDS, LIQUIDS OR PHARMACEUTICAL PRODUCTS.

Important notes



THE INFORMATION INCLUDED IN THIS MANUAL ARE SUBJECT TO CHANGES WITHOUT ANY NOTICE.

MO.COM. LTD. CO. WON'T BE RESPONSIBLE FOR DIRECT, INDIRECT, ACCIDENTAL, CONSEQUENT DAMAGES OR OTHER DAMAGES RELATED TO THE SUPPLY OR THE USE OF SUCH INFORMATION.

ANY REPRODUCTION, MODIFICATIONS OR TRANSLATION OF THIS MANUAL OR PARTS OF IT IS FORBIDDEN WITHOUT PREVIOUS WRITTEN AUTHORIZATION BY M.O.COM. LTD. CO.

mocom AND **millennium** ARE REGISTERED TRADEMARKS OF M.O.COM. SRL.

SCOPE OF THE MANUAL

The present manual has the scope to supply instructions for:

- general knowledge of the product;
- correct installation and configuration;
- sure and efficient use;
- correct treatment of the material before and after the sterilization.

In the appendix you will find:

- technical characteristics of the product;
- specifications of the sterilization PROGRAMMES;
- maintenance procedures;
- analysis and solution of the problems;
- additional documentation.

GENERAL WARNINGS

The product should be used in compliance with the procedures described in the manual and never for purposes differing from the foreseen ones.



THE USER IS RESPONSIBLE FOR THE LEGAL FULFILLMENTS CONCERNING THE INSTALLATION AND USE OF THE PRODUCT.

IF THE PRODUCT IS NOT CORRECTLY INSTALLED AND USED, OR A SUITABLE MAINTENANCE IS NOT OPERATED, THE MANUFACTURER CANNOT BE CONSIDERED RESPONSIBLE OF POSSIBLE BREAKS, MALFUNCTIONS, DAMAGES, LESIONS TO PROPERTY AND/OR PEOPLE.

In order to avoid any danger situations, with possible consequent damages or lesions to property and/or people, the following precautions should be followed:

- Use high quality distilled water only.





THE USE OF WATER OF INADEQUATE QUALITY COULD SERIOUSLY DAMAGE THE EQUIPMENT.

SEE APPENDIX A, TECHNICAL CHARACTERISTICS.

- Do not pour any water or liquids over the equipment;
- Do not pour inflammable substances over the equipment;
- Do not operate on the equipment in presence of explosive or inflammable gas;
- Before any maintenance or cleaning action **ALWAYS REMOVE** the power mains.

If this precaution is impossible or the external mains breaker is far or not visible from the people performing the maintenance, affix the poster **WORKS IN PROGRESS** on the external breaker after having positioned it on OFF;



- **Make sure** that the electric plant is provided with the earth connection in accordance with the current laws;
- **Do not remove** any label or plate from the equipment; in case call for new ones.
- **Use exclusively** original spare parts.




THE NOT OBSERVANCE OF WHAT ABOVE DESCRIBED MAKES ANY RESPONSIBILITY OF THE MANUFACTURER TO DECAY .

PACKING CONTENT

DIMENSIONS AND WEIGHT

Height	490 mm
Width	800 mm
Depth	485 mm

Total weight ~ 50 kg



ON RECEIVING THE PRODUCT, PLEASE VERIFY THAT THE PACKING IS INTACT IN EVERY PART.


Open the pack and check that:

- The supply corresponds to the order specifications (see delivery note);
- there are not evident damages on the product.



IN CASE OF WRONG SUPPLY, MISSING PARTS OR DAMAGES OF WHATEVER TYPE, PLEASE IMMEDIATELY INFORM AND IN DETAIL THE RETAILER AND THE FORWARDER THAT CARRIED OUT THE DELIVERY.

CONTENT



The packing contains the following items:

- **millennium B_μ** steam sterilizer;
- Two stainless steel wire trays (Ref. 1);
- One stand stainless steel tray-holders (Ref. 2);
- Instrument-holder box (Ref. 3);
- Stainless steel box holder (Ref. 4)
- Tray removal tong (Ref. 5);
- 2-liters capacity tank provided with fast fitting-ended hose for manual filling (Ref. 6);
- Additional bacteriological filter (Ref. 7).
- Draining silicon hose (1 m), ended with fast fitting (Ref. 8).
- Operating manual (Ref. 9)
- Warranty certificate (Ref. 10) (see note).



KEEP HOME THE GUARANTEE CERTIFICATE TOGETHER THE INVOICE.

**LIFTING AND
TRANSPORT**

The packed product should be handled by using, where possible, suitable mechanic tools (cart elevator, transpallet, etc.) and following the indications printed on the packing.
In case of manual handling, the product should be moved by two people and using the proper handles provided on the box.

The sterilizer should be lifted out the box by two people and moved through a cart or similar means.



WE RECOMMEND TO TRANSPORT AND STORE THE EQUIPMENT AT A TEMPERATURE HIGHER THAN 5°C. A PROLONGED EXPOSURE AT LOWER TEMPERATURE COULD CAUSE DAMAGES TO THE PRODUCT .



STORE THE ORIGINAL PACK THAT WILL BE USED FOR A POSSIBLE FUTURE TRANSPORT OF THE EQUIPMENT. USE OF DIFFERENT PACK COULD CAUSE DAMAGES TO THE PRODUCT ON SHIPPING.



BEFORE MOVING THE STERILIZER IT IS NECESSARY TO EMPTY THE DISTILLED AND RECOVERY WATER TANKS, AFTER THE EQUIPMENT HAS BEEN SWITCHED OFF FOR ABOUT 60 MINUTES FROM THE LAST CYCLE, IN ORDER TO ALLOW THE COOLING OF ALL THE INTERNAL HOT ELEMENTS.

GENERAL DESCRIPTION

INTRODUCTION

Millennium B_μ is the revolutionary MO.COM. proposal and represents the "state of art" concerning speed, safety and performances as well as the technological frontier in the field of the fast and small B-type (prEN 13060: 2004) steam sterilizers.

This sterilizer is a sophisticated equipment but of very friendly use, adaptable to the different demands thanks to the wide possibilities of configuration and choice of the cycles.

The equipment is able to treat in fast way every type of load, in particular handpieces, turbines and hollow articles, through various patented systems and the operation completely microprocessor controlled.

Besides, this equipment allows a better approach by the user, that, rather than to conform himself to the machine characteristics, can "communicate" with it for configuring the performances to the different job requirements.

Thanks to the friendly use, dimensions very reduced and unique agreeable design, it represents the ideal partner for the professional requiring the maximum sterilization safety.

GENERAL FEATURES

Millennium B_μ is a microprocessor-based water steam sterilizer, equipped with a new-type ovalized sterilization chamber made of anticorrosive extrusion (volume of 5.5 litres).

It is characterized by an advanced fractionated vacuum system for the complete air removal also from hollow and porous materials as well as by an effective final vacuum drying phase able to eliminate any trace of condensation from the load.

An exclusive steam generation system, combined with the electronic management and high accuracy sensors, guarantees an high process speed and high stability of the thermodynamic parameters during the whole sterilization process.

Process Evaluation System besides monitors in real time all "key" parameters of the equipment in order to warranty full safety and perfect result.

The equipment offers 11 sterilization PROGRAMMES (including one program completely programmable), optimized for an effective and fast sterilization of different tools and materials used in medical environment, particularly the dental one.

Four PROGRAMMES can be directly recalled via the command panel, characterized with simple and new agreeable layout design.



Besides, the sterilizer offers the facility to choice the preheating modes according to the usage frequency, the print options of the cycle report, the water filling and draining modes and other features.

For more information, refer to Chapter "**Setting**".


Finally **Millennium B_μ** has the most sophisticated and advanced safety systems in order to guarantee the user against possible operation anomalies, both electric, mechanics, thermal and biologic type.




FOR THE DESCRIPTION OF THE SAFETY DEVICES REFER TO APPENDIX A (TECHNICAL CHARACTERISTICS).

FRONT DEVICES**REAR DEVICES**

COMMAND PANEL



LCD DISPLAY



**OPERATING
CYCLE EXAMPLE**

The sterilization program can properly be described by a sequence of phases, each with a well defined activity.


Considering the example of the standard program (i.e. program at 134°C - 4' for porous material), after arranged the material in the chamber, closed the door, selected the program and started the cycle by the START key, and consequent activation of the door locking mechanism, the following sequence will run:

1. Pre-heating of the steam generator and sterilization chamber;
2. Chamber air removal and steam penetration inside the load through a series of vacuum (fluid exhaustion from the sterilization chamber) and pressure phases (steam entering into the chamber);
3. Pressure increasing, with consequent steam temperature increasing up to the preset sterilization conditions (for the example, **134°C**);
4. Stabilization of the pressure and temperature conditions inside the sterilization chamber;
5. Run of the sterilization process for the preset time (for the example, **4 minutes**);
6. Chamber pressure decreasing through steam discharge;
7. Vacuum drying phase;
8. Venting phase through sterile air;
9. Chamber pressure levelling up to atmospheric value.

As reached this last phase the door locking mechanism will be released and the door can be opened to recovery the load from the sterilization chamber paying attention to use proper precautions.

Whereas the phases 1, 3, 4, 6 and 9 are essentially identical for all the cycle types, with little time differences depending on the quantity and consistence of the load and heating conditions, the phases 2, 5, 7 and 8 vary their configuration and/or duration sharply depending on the selected cycle (and accordingly on the load typology) and on the choices operated by the user.

A standard program can be outlined as the following graph.



FOR DETAILS ABOUT THE AVAILABLE PROGRAMMES, REFER TO APPENDIX **B** (PROGRAMMES).

INSTALLATION

GENERAL

For the right operation, features exploitation and life extension of the sterilizer, the first and basic step is the correct and careful installation. Such precaution besides avoids possible malfunctions or damages to the equipment, as well as possible danger situations for property and people.


We suggest to **meticulously** follow the warnings reported on this chapter.

☞ “ASSISTENZA CLIENTI M.O.COM.” (SEE APPENDIX Z) IS AT YOUR DISPOSAL FOR ANY QUESTION OR FURTHER INFORMATION.

☞ THE STERILIZER IS FORWARDED AFTER HAVING PASSED A SET OF PROGRAMMED CONTROLS. THEREFORE FURTHER CALIBRATIONS FOR THE OPERATION ARE NOT NECESSARY ANY MORE

Dimensions and weight

– Height (total)	340 mm
– Width (total)	380 mm
– Total depth (excluded rear fittings)	585 mm
– Depth for build mounting	470 mm
– Weight (total)	38 kg



Power supply

The electric plant to which the sterilizer is connected has to correctly be sized according to the electric characteristics of the equipment. The rating data are reported on the plate on the **rear panel** of the machine.

☞ THE MAINS SOCKET SHOULD BE EASILY ACCESSIBLE AND NEAR THE EQUIPMENT.



VERIFY THE ADEQUACY OF THE ELECTRIC PLANT, PARTICULARLY FOR THE EARTH CONNECTION.

DIMENSIONS FOR BUILD MOUNTING

In case the sterilizer has to be set into a piece of furniture it is necessary to provide a suitable space all around the equipment in order to assure an effective ventilation, as well as a free space in the rear in order to allow the passage of the draining pipelines and power supply cord, and to assure an adequate air flow for the cooling of the heating exchanger.

☞ LEAVE A SUITABLE SPACE FOR MAINTENING PURPOSE OF THE BACTERIOLOGIC FILTER LOCATED ON THE REAR OF THE EQUIPMENT. PERIODICALLY PERFORM THE CHECK OF THE FILTER WEAR AND TEAR (SEE APPENDIX C, MAINTENANCE).

Accordingly it is important that the build for the installation has the **minimum dimensions** as in figure:



BUILD DIMENSIONS SMALLER THAN THE INDICATED ONES CAN JEOPARDIZE THE AIR CIRCULATION AROUND THE EQUIPMENT AND NOT GUARANTEE A SUITABLE COOLING, WITH CONSEQUENT REDUCTION OF THE PERFORMANCES AND/OR POSSIBLE DAMAGES.




IN CASE THE STERILIZER INSTALLATION MAKES THE MAINS BREAKER INACCESSIBLE, USE A MAINS SOCKET INCORPORATING A SWITCH.



DON'T REMOVE THE CARTER NEITHER OTHER EXTERNAL ELEMENTS. ARRANGE INTO THE BUILD THE EQUIPMENT COMPLETE IN EVERY PART.



FOR THE COMPLETE TECHNICAL DATA REFER TO APPENDIX A. (TECHNICAL CHARACTERISTICS).



GENERAL INSTALLATION DIRECTIONS

In order to assure the correct operation of the equipment and/or avoid situations of risk, observe the **following directions**:

- Install the sterilizer on a plain surface; eventually adjust the rear feet to match possible irregularities.
Verify that the plan is able to support the weight of the equipment (~ 45 kg);
- **Allow a suitable space for the ventilation (at least 10 cm each side) all around the sterilizer, particularly on the back side.**
In case the equipment is set into a piece of furniture, verify the respect of the warnings as above indicated, avoiding possible obstruction of the air slots;
- Do not install the sterilizer near tubs, sinks or analogous places, avoiding the contact with water or liquid. Otherwise, short circuits and/or situations of potential danger for the operator could occur;
- Do not install the sterilizer in environments characterized by the presence of excessive humidity or in rooms with poor ventilation;
- Do not install the sterilizer in environments characterized by the presence of inflammable and explosive gas or vapours;
- Install the equipment in such a way the power supply cable could not be bent or crushed. The cable should run freely to A.C. socket.



INSTALL THE EQUIPMENT IN SUCH A WAY THE POSSIBLE EXTERNAL LOAD/DRAIN PIPELINES COULD NOT BE BENT OR CRUSHED. THEY SHOULD RUN FREELY TO THE EXTERNAL DRAINING TANK.


MAINS CONNECTION

The sterilizer should be connected, in compliance with the laws and/or directives in-force, to a mains socket of the electric plant with a power rate suitable for the equipment consumption, provided with earth connection and max impedance of: |Z line| = 0.0109; |Z neutral| = 0.0073.

The mains socket should be suitably protected through differential switch having the following characteristics:

- Rating current I_n **16 A**
- Differential current $I_{\Delta n}$ **0,03 A**

Differential switch



Mains socket with earth pin



Mains switch



THE MANUFACTURER IS NOT LIABLE FOR THE DAMAGES CAUSED BY AN INSTALLATION WHERE A NOT SUITABLE ELECTRIC PLANT AND/OR EARTH CONNECTION ARE PROVIDED.

In case the mains plug doesn't match the socket, replace the cable plug with a suitable type of same electrical characteristics or anyway suitable for the electrical requirements of the equipment. The plug choice and replacement is on care and responsibility of the user.




DIRECTLY CONNECT THE CABLE TO THE MAINS SOCKET. **DON'T USE** EXTENSIONS, ADAPTERS OR OTHER ACCESSORIES.

EXTERNAL PRINTER CONNECTION (OPTION)

Connect the printer parallel cable to DB-25 female connector on rear side of the sterilizer.

Any type of printer can be supported (impact, thermal, ink-jet or laser printer); in any case verify the correct operation of the sterilizer with the printer selected.

For the printer parameter setting, refer to chapter "**Setting the equipment – Setting the printer option**".



FIRST START-UP

Once the sterilizer has correctly been installed, you can proceed to the switching on and equipment configuration.

SWITCHING ON

Switch on the sterilizer through the mains switch mounted on the rear.




PERFORM THIS ACTION WITH THE DOOR OF THE STERILIZER IN OPEN POSITION.

INITIAL SELF-TEST

At the switching on, the equipment turns on the LCD display and generates a signaling tone, allowing the user a visual check of the correct command panel operation.


Then the display will show:



IN CASE OF CLOSED DOOR, THE SEL-TEST IS STOPPED, A WARNING TONE IS GENERATED AND THE FOLLOWING MESSAGE ON THE LCD DISPLAY:




Open the door to allow the test continuing. Over the test the following message is shown:



ACQUISITION AND AUTOMATIC UPDATING OF THE ENVIRONMENT PRESSURE VALUE

For the correct working of internal auxiliary devices the unit, at the switching on, provides for the automatic reading of the environmental pressure value. If the difference between the currently read pressure value and the pressure value stored (see Chapter “**Setting the equipment**”, par. **Setting the options in SETUP mode**) is higher than a preset value, the unit will update automatically the stored value after a short waiting time, otherwise the stored value remains unchanged.

After this check and in case of updating, the unit starts the self-test (see previous paragraph), and the display shows the following message with a warning tone at the end:




Press the key ↓ to enter in STAND-BY status (see the following).



FOR MORE INFORMATION SEE THE PARAGRAPH “ACQUISITION OF THE ENVIRONMENT PRESSURE VALUE”.

STAND-BY STATUS

Over the self-test procedure the sterilizer enters in **STAND-BY** mode, and the LCD display will show:



where the first line is showing the number of the cycles correctly performed (xxxxx) and the number of the total cycles launched (yyyyy), second line the Stand-by status (High-Low-Off), next line the current temperature value and **date**, and the last line the current pressure values and **time**.




THE CYCLE IS SEEN AS LAUNCHED WHEN STARTS THE FIRST VACUUM PHASE, THE PREHEATING PHASE EXCLUDED. THE CYCLE IS SEEN AS COMPLETED WHEN REACHES THE END OF THE PROGRAM (SEE CHAPTER "RUNNING THE STERILIZATION PROGRAM").



FOR DATE AND TIME SETTING AS WELL AS THE CHOICE OF THE PREHEATING, PRINT AND FILLING MODES REFER TO CHAPTER "SETTING THE EQUIPMENT".

At regular interval, the two upper lines will show alternatively the above content and the current setting for print (Off-On) and filling (Manual) mode.



The icons on the lower part of the LCD display are off, except the door status and tank level indications that will be on respectively if the door is closed or in case of MIN/MAX water level detection in the distilled water tank (or MAX water level detection in the recovery tank).

At the first start-up, the icon relating the MIN water level of the distilled water reservoir will be normally on.

The equipment is waiting for the selection of the sterilization program by the user (see Chapter "Selecting the sterilization program").




HOT SURFACES: IN STAND-BY MODE AND THE DOOR OPEN, A BEEP SIGNAL WARNS FOR HIGH TEMPERATURE OF THE INTERNAL SURFACES. IN ORDER TO AVOID SCALDS, DO NOT TOUCH BY HANDS THE STERILIZATION CHAMBER, THE TRAYS OR THE INNER PART OF THE DOOR.

FILLING THE DISTILLED WATER

Before using the sterilizer at the first start-up or any time the MIN water level signaling turns on, it is necessary to fill or top-up the distilled water reservoir.

Operate in the following way (with door open) and refer to the figure:

1. Hold horizontally the manual filling tank and fill it of distilled water (2 l);
2. Connect the fast fitting of the supplied pipe to the inlet mounted on the front upper side of the equipment, pushing for a click;
3. Position vertically the tank and loosen the cap taking care do not upset water on the machine.
4. Now the water will start to flow into the internal reservoir;
5. On continuing the filling, the MIN level indicator turns off.
6. Continue up to dry out completely the tank;
7. Take the tank and lower it below the connection point, holding it horizontally;
8. Push on the metallic clip of the fitting and remove the rubber pipe;
9. Fill again the tank (2 l) and repeat once again the steps 2, 3 and 4;
10. As the icon **MAX** signaling turns on (accompanied by an acoustic signal) interrupt the filling and operates as described at the steps 7 and 8.



TO START THE STERILIZATION PROGRAM IS NOT NECESSARY THAT THE **MAX** SIGNALING IS TURNED ON. IT IS SUFFICIENT THAT THE **MIN** SIGNALING IS TURNED OFF




USE ONLY HIGH QUALITY DISTILLED WATER. FOR THE INDICATIONS CONCERNING THE WATER FOR STERILIZATION EQUIPMENT SEE APPENDIX A (TECHNICAL CHARACTERISTICS).

RECOVERY TANK MAX WATER LEVEL SIGNALING

The maximum water level in the recovery tank is properly signaled by the blinking of the MIN and MAX indicator on LCD display.

In this case provide for the draining of the internal recovery tank as follows:

- Arrange an external tank with capacity of at least 4 liters;
- Open the door and position the free side of the supplied tube in the external tank and the plug the opposite side provided with fast coupling fitting in the draining point mounted on the front bottom side of the equipment, pushing for a click;
- Drain completely the water recovery reservoir;
- Push on the metallic clip of the fitting and remove the rubber pipe.



SETTING THE EQUIPMENT

INTRODUCTION

millennium B_μ allows new and wider customization facilities never offered before by any sterilizer of its category.

The configuration of the sterilizer can be set according with the user's requirements, by adapting, for instance, the characteristics of the equipment depending on the working frequency and the type of material to be treated.

Through a sophisticated SETUP program a lot of options can be selected by the user through a friendly and easy menu.




USE THE **SETUP** PROGRAM EVERY TIME YOU NEED. A CORRECT CUSTOMIZATION OF THE EQUIPMENT ALLOWS TO GET THE BEST PERFORMANCES AND THE MAXIMUM SATISFACTION FOR THE USE.




"ASSISTENZA CLIENTI M.O.COM." (SEE APPENDIX Z) IS AT DISPOSAL OF THE USERS TO SUPPLY SUGGESTIONS OR ADVICES FOR THE BETTER USE THE OPTIONS AVAILABLE IN THE **SETUP** PROGRAM.

STARTING THE SETUP PROGRAM



To enable the **SETUP** program push for a few seconds on the **SETUP** key (↑) on the command panel, up to the display is showing:



ICON **SETUP** ON LCD TURNS ON FOR THE TIME THE SETUP MODE IS ENABLED.

Press the key ↓ to enter the SETUP program. Now the display is showing the first-level menu options (see the **Flowchart of the SETUP program**).

On the contrary, press the key **ESC** ↑ to exit the SETUP program and return in the stand-by mode.




THE **SETUP** PROGRAM CAN ONLY BE RECALLED FROM **STAND-BY** MODE. DURING THE STERILIZATION PROGRAM THE **SETUP** MODE CANNOT BE ACCESSED.

SETUP MODE KEYS FUNCTION

Differing from the normal operation, in **SETUP** mode the keys have the function corresponding to symbol marked over the key as follows:

- Symbol ↓ **ENTER** function, to confirm the data (START/STOP key)
- Symbol + **increment/scrolling up** function (program selection key)
- Symbol — **decrement/scrolling down** function (test selection key)
- Symbol ↑ **ESC** function, to exit the option (SETUP key)

These function keys are valid as long as you remain in the **SETUP** mode.



DESCRIPTION OF THE MENU OPTIONS

On the following the meaning of the main menu and second-level options are described.

MAIN menu

The main menu of the SETUP program has 6 options that recall further menus (second-level):

BASIC	(basic options)
ADVANCED	(advanced options)
SPECIAL	(special options)
SERVICE	(<i>not user-available</i>)
DATA REVIEW	(review of the selected options)
EXIT SETUP	(exit the SETUP mode and return to normal operation. See paragraph Exit the SETUP mode)



THE WAY FOR SETTING THE DIFFERENT OPTIONS IS DESCRIBED ON THE PARAGRAPH "**SETTING THE OPTIONS IN THE SETUP MODE**".

BASIC menu

This menu consists of the following options:

LANGUAGE	(current language setting)
DATE SETTING	(current date setting);
TIME SETTING	(current time setting)
PASSWORD	(password setting)
EXIT	(exit the BASIC menu and return to main menu)

ADVANCED menu

It consists of the following options:

PROGRAMMES	(preset sterilization programmes selection, as shown on LCD)
STAND-BY OPTIONS	(stand-by mode setting)
PRINT OPTIONS	(printer and printing option setting)
FILLING OPTIONS	(not applicable)
DRAIN OPTIONS	(not applicable)
EXIT	(exit the ADVANCED menu and return to main menu)

SPECIAL menu

This menu consists of the following options:

AMBIENT PRESSURE	(ambient pressure acquiring)
LCD CONTRAST	(display contrast adjustment)
EXIT	(exit the SPECIAL menu and return to main menu)

SERVICE menu

This menu can be accessed by Customer Service **ONLY**.

DATA REVIEW menu

This menu allows to recall the **current settings** of the equipment, allowing the user to verify the exactness.

The following screens are available (example of current setup).

M I L L E N N I U M B μ R . E x x x x / M M y y y y y y L A N G U A G E E N G L I S H					D A T E d d / m m / y y y y T I M E h h : m m : s s				
			SETUP					SETUP	

Firmware version

Use the keys + / - to scroll.

1 s t P R E S E T 1 3 4 ° C P O R O U S S T A N D A R D D R Y I N G					2 n d P R E S E T 1 3 4 ° C H O L L O W F A S T D R Y I N G				
			SETUP					SETUP	

Use the keys + / - to scroll.

3 r d P R E S E T 1 3 4 ° C W R A P P E D E X T R A D R Y I N G + 0 5					4 t h P R E S E T 1 3 4 ° C E M E R G E N C Y F A S T D R Y I N G				
			SETUP					SETUP	

Use the keys + / - to scroll.

S t a n d - b y o p t i o n H I G H 1 2 0 m i n					P r i n t o p t i o n O F F 1 C O P Y (i e s)				
			SETUP					SETUP	

Use the keys + / - to scroll.

F i l l i n g o p t i o n M A N U A L D r a i n o p t i o n I N T E R N A L				
			SETUP	

Use the keys + / - to scroll.

E X I T D A T A R E V I E W ↓ t o c o n t i n u e				
			SETUP	

Press ↓ to confirm.



FOR THE MEANING OF THE ABOVE TERMS SEE PARAGRAPH "SETTING THE OPTIONS IN SETUP MODE".

DEFAULT SETTINGS

millennium B_μ is delivered with the following default settings:

DATE: *current date*
TIME: *current time*

PROGRAMMES: Preset 1: **134°C POROUS** (*standard* drying)
 Preset 2: **134°C HOLLOW** (*standard* drying)
 Preset 3: **134°C SOLID** (*standard* drying)
 Preset 4: **134°C EMERGENCY**



THE ABOVE PROGRAMMES ARE PREFERRED SETTINGS. ANY COMBINATIONS ARE POSSIBLE FOR DIFFERENT MARKETS.

ST-BY MODE: **HIGH** (pre-heating)

PRINT OPTIONS: **OFF**

FILLING OPTIONS: **MANUAL**

DRAIN OPTIONS: **INTERNAL**

SETTING THE OPTIONS IN SETUP MODE

Language

(LANGUAGE option in BASIC menu)

We analyze now in detail how to select the different options, proceeding with the order as shown in the previous paragraph (Description of the menu options).

Select the option **LANGUAGE** with the key \downarrow , the following screen will be shown:



Select the desired language.

Scroll the list by the key $+$ or $-$ and confirm with key \downarrow to store the choice.

As confirmed the choice you re-enter in the second-level menu.




FROM NOW ON, THE SETUP MENU WILL BE SHOWN IN THE LANGUAGE SELECTED.

Date setting

(DATE SETTING option in BASIC menu)

As selected the option **DATE SETTING** by key \downarrow , the following screen will be shown:



Carry out the following:

- Day indication **is flashing**: adjust for the current date by keys $+$ and $-$. Confirm by \downarrow .
- Month indication **is flashing**: adjust for the current month by keys $+$ and $-$. Confirm by \downarrow .
- Year indication **is flashing**: adjust for the current year by the keys $+$ and $-$. Confirm by \downarrow .

As entered the last item the date will be stored and the second level of the menu returned.

Time setting

(TIME SETTING option in BASIC menu)

Selecting the option **TIME SETTING** by key ↵, the following screen will be shown:

h h : m m : s s									
+	/	-	t	o	s	e	t		
↵	t	o	e	n	t	e	r		
↑	t	o	e	x	i	t			
							SETUP		

Carry out the following:

- Hour indication **is flashing**: adjust for the current hour by keys + and -. Confirm by ↵.
- Minute indication **is flashing**: adjust for the current minutes by keys + and -. Confirm by ↵.

As entered the last item the time will be stored and the second level of the menu returned.

Password setting

(PASSWORD option in BASIC menu)

Selecting the option **PASSWORD** by key ↵, the next screen will be shown:

→	D I S A B L E D										+
	A	N	Y	P	O	W	E	R	O	N	↑
	A	N	Y	C	Y	C	L	E	S	T	↓
	E	X	I	T							-
							SETUP				

Select **DISABLED** to use the sterilizer freely without no restriction for the operator access.

Select **ANY POWER ON** to safeguard by password the sterilizer from undesired switching on (through mains switch).

In this way the sterilizer will be switched on only by authorized personnel and then used freely by any people.

Select **ANY CYCLE START** to safeguard by password the sterilizer from undesired switching on and use.

Only authorized personnel will be allowed to power on the equipment and enter a command to start a sterilization cycle.



THIS FEATURE WARRANTS A MORE CONTROLLED USE OF THE STERILIZER, EVEN IF THE PASSWORD MANAGEMENT RESULTS MORE COMPLICATED. ENABLE THIS FEATURE ONLY IN CASE OF ACTUAL NEED TO AVOID UNNECESSARY COMPLICATION.

Selecting an enabling password option a next screen will be shown:

I N S E R T P A S S W O R D											
↵	t	o	e	n	t	e	r				
↑	t	o	e	x	i	t					
							SETUP				

Enter the password (8 characters, fixed length) by the keys + and - and confirm by the key ↵. The same screen is shown again:

C O N F I R M P A S S W O R D											
↵	t	o	e	n	t	e	r				
↑	t	o	e	x	i	t					
							SETUP				

Reenter the password and confirm by the key ↵.



TO MODIFY THE PASSWORD, FIRSTLY SELECT THE OPTION DISABLED TO DELETE THE PASSWORD PREVIOUSLY ENTERED, THEN SELECT THE OPTION ANY POWER ON OR ANY CYCLE START AND ENTER THE NEW PASSWORD AS ABOVE INDICATED.

Setting the preset sterilization Programmes

(PROGRAMMES option in ADVANCED menu)

The selection of the programmes and their storing under a preset position should be performed by the following steps through different menu sequence.

It is possible to choose both the **preset** PROGRAMMES and the user **configurable** program (CUSTOM).

We will see the two cases.

To associate a **PRESET** program to a position operate as follows:

- As selected the **PROGRAMMES** option by key ↵, the following menu appears:

→	1 s t	P R E S E T	+
	2 n d	P R E S E T	↑
	3 r d	P R E S E T	↓
	4 t h	P R E S E T	-
	E X I T		

Define the position (1, 2, 3 or 4) the preset sterilization program has to be associated and select the item by keys + and -. Confirm by ↵.

- Now the display is showing the list of the available cycles (four lines in turn):

→	1 3 4 ° C	H O L L O W	+
	1 2 1 ° C	H O L L O W	↑
	1 3 4 ° C	S O L I D	↓
	1 2 1 ° C	S O L I D	-
.....			
	1 3 4 ° C	E M E R G E N C Y	
	1 3 4 ° C	W R A P P E D	
	1 2 1 ° C	W R A P P E D	
	1 3 4 ° C	P O R O U S	
	1 2 1 ° C	P R O R O U S	
	1 3 4 ° C	P R I O N	
	X X X ° C	C U S T O M	
	E X I T		

Scroll the list through the keys + and - to find the desired sterilization program.

- Enter the choice by ↵.

In case of **PRION** choice, the following screen is proposed to select the sterilization time.

You can set a time value starting from **18** minutes.

Depending on the choice performed, two alternative menus are offered to select the drying mode to be associated with the program selected .

a) PROGRAMMES with short drying (HOLLOW, SOLID, EMERGENCY):

You can choice the mode **STANDARD** (default) or **FAST** (short drying, preferred with light load).

Move the cursor by key + and - . Confirm by ↵.



EMERGENCY PROGRAM ALLOWS ONLY THE FAST DRYING OPTION.

b) PROGRAMMES with long drying (POROUS, WRAPPED, PRION):

You can select **STANDARD** (default), **INTELLIGENT** (automatic drying characterized with a duration - higher or lower the standard one - depending on the volume, quantity and typology of the load) or **EXTRA** (drying phase extended of a selectable time period - suggested for particularly difficult load

Use keys + / - and confirm the choice by key ↵.



WITH BIG LOAD OR SPECIAL MATERIAL, THE **STANDARD** OPTION MIGHT **NOT** SOMETIMES GUARANTEE A CORRECT RESULT. IN THIS CASE, EXTEND THE DRYING PHASE WITH THE **EXTRA** OPTION.



WITH COMPLEX LOAD (E.G. WRAPPED INSTRUMENTS IN A STERILIZATION “CONTAINER”) THE “INTELLIGENT” OPTION MIGHT NOT CORRECTLY OPERATE, WITH RESULTS LOWER THE EXPECTATION. IN THIS CASE USE THE **STANDARD OR EXTRA** OPTIONS AS NECESSARY.

Enabling the **EXTRA** option, the following choice will be prompted:

allowing to set the time period of the extra drying.

The extended time ranges from **1** to **4** minutes (added to the STANDARD drying time)


Use keys + / - to set the time and confirm the choice by key ↵.



THE SELECTION CAN BE WHENEVER MODIFIABLE BY PERFORMING AGAIN THE PROCEDURE ABOVE DESCRIBED.




IF AN IDENTICAL STERILIZATION PROGRAM IS ALREADY PRESENT IN A DIFFERENT PRESET POSITION, THE CHOICE WON'T BE ACCEPTED. ON THE SCREEN THE FOLLOWING WARNING APPEARS, ACCOMPANIED BY ACOUSTIC SIGNALINGS:




To configure a **CUSTOM** cycle and associate it to a position (1, 2, 3 or 4), proceed as follows:

1. Select the item **PROGRAMMES**, select the number to which the program has to be associated (as previously described), then select the item **CUSTOM** in the next screen; the following menu appears:



Select **121°C** for a custom sterilization program at 121°C or **134°C** for a program at 134°C. Move by keys + and -. Confirm by ↵.

2. The next screen:




requiring to choose the duration of the sterilization process. Use the key + / - to adjust the value, then confirm by the key ↵.



THE TIME OF THE STERILIZATION PROCESS RANGES BETWEEN 4 AND 30 MINUTES FOR THE PROGRAM AT 134°C, AND BETWEEN 20 AND 30 MINUTES FOR THE PROGRAM AT 121°C.


3. Next, appears the menu for selecting the type of the initial vacuum:



Select **FRACTION. VACUUM** to perform a fractionated vacuum (suggested for hollow and porous materials), or **SINGLE VACUUM** for one pre-vacuum phase only (solid instruments).

Use the key + / - then confirm by the key ↵.


4. Now a next menu is shown where the drying modes can be selected:



Select **LONG DRYING** to properly dry porous and/or wrapped instruments or **SHORT DRYING** in case of solid (or hollow) unwrapped instruments.
Use key + / - then confirm by key ↵.

5. Depending on the option selected (**SHORT** or **LONG**), two different menus will be prompted.

LONG DRYING:



Select **STANDARD DRYING** for a drying with fixed time.


- ☞ WITH BIG LOAD OR SPECIAL MATERIAL, THE **STANDARD** OPTION MIGHT **NOT** GUARANTEE SOMETIMES A CORRECT RESULT. IN THIS CASE, EXTEND THE DRYING PHASE WITH THE **EXTRA** OPTION.

Select **INTELL. DRYING** to enable an automatic drying characterized by a duration (higher or lower the standard one) depending on load volume, quantity and typology.

- ☞ WITH COMPLEX LOAD (E.G. WRAPPED INSTRUMENTS IN A STERILIZATION “CONTAINER”) THE “INTELLIGENT” OPTION MIGHT NOT CORRECTLY OPERATE, WITH RESULTS LOWER THE EXPECTATION. IN THIS CASE USE **STANDARD** OR **EXTRA** OPTIONS AS NECESSARY.

Select **EXTRA** option for extending the drying phase of a selectable time period (e.g. for particularly difficult load).

Enabling this option, the following choice will be prompted:




allowing to set the time period of the extra drying.

The extended time ranges from **1** to **4** minutes (added to the **STANDARD** drying time)

Use the keys + / - to set the time and confirm the choice by key ↵.

- ☞ IF THE **CUSTOM** PROGRAM IS ALREADY PRESENT IN A DIFFERENT PRESET POSITION, THE CHOICE WON'T BE ACCEPTED. ON THE SCREEN THE FOLLOWING WARNING APPEARS, ACCOMPANIED BY ACOUSTIC SIGNALING:



- ☞ THE SELECTION CAN BE WHENEVER MODIFIABLE BY PERFORMING AGAIN THE PROCEDURE ABOVE DESCRIBED.

- ☞ THE LIST OF THE AVAILABLE PROGRAMMES, THEIR DIAGRAMS AND THE CHARACTERISTICS OF THE AUTOCLAVABLE MATERIAL (DEPENDING ON THE PROGRAM) CAN BE FOUND IN **APPENDIX B (PROGRAMMES)**.


- ☞ NO PASSWORD IS REQUIRED TO ACCESS THE **USER** PROGRAMME. ANY POSSIBLE COMBINATIONS OF THE CYCLE PARAMETER CUSTOMIZATION DO NOT INVOLVE HAZARDS TO USER OR EQUIPMENT.

Setting the STAND-BY modes


(STAND-BY OPTIONS in the ADVANCED menu)

Value and timeout of the pre-heating phase in STAND-BY mode can be set according to the high- or low-use of the equipment and others possible considerations.

Selecting **STAND-BY OPTIONS** by key ↵, the following menu will appears:



Selecting **ST-BY MODE** option a next menu is offered for setting the heating level:



Select **HIGH** (high preheating) for intensive use and to reduce however the waiting time between one cycle and the other.

Select **LOW** (low preheating) for normal use, maintaining an acceptable waiting time.


Select **OFF** (preheating switched off) for non-continuous use.

In this case the waiting time will be longer (8-10 minutes starting from the “cold” start condition).

Use the key + / - then confirm by the key ↵.

Selecting **ST-BY TIME-OUT** option, you can set the timeout for disabling the STAND-BY mode, i.e. the time after the last cycle the heating system must be disabled.

The following screen appears:



You can range from **0** to **300** minutes (in step of 30 minutes), after which the heater will be disabled (as stand-by mode in OFF condition) to avoid useless power consumption.

Adjust the value by means the keys + / - then confirm with the key ↵.



THIS FUNCTION IS ALSO ACTIVE WITH STAND-BY SET TO OFF, HOWEVER THE TIMEOUT HAS NO MEANING AS THE HEATER IS TURNED OFF AT THE END OF THE CYCLE..



PRESSING ANY CYCLE KEY (STERILIZATION OR TEST KEY), OR SWITCHING OFF/ON THE UNIT BY THE MAINS SWITCH, THE ORIGINAL STAND-BY MODE (HIGH OR LOW) WILL BE RESTORED.

Setting the Printer Options

(PRINT OPTIONS of the ADVANCED menu)

In case the sterilizer is connected to a printer for the recording of the sterilization data, set the necessary parameters in order to guarantee the correct operation.

Proceed as follows:

1. By selecting the item **PRINT OPTIONS** and confirming it by the key ↵, appear the following menu:



The item **PRINTER** opens a next menu for selecting the line feed mode accordingly with the printer used, whereas the item **REPORT** is used to set the number of the printouts and to enable the reprint of the last program performed.

a) PRINTER option


The following screen is proposed:



Choice **OFF** to disable the printing after the cycle (sterilization or test cycle).

In **millennium B_μ** the option **INTERNAL** is not used. If selected, the sterilizer will continue with the option **EXTERNAL**.

Select **EXTERNAL** to enable the printing of the data on an external peripheral. The following menu appears:



Choice **CR** if you are using printers that carry out the form feed through a CR command (*Carriage Return*), or choice **CR+LF** for printers requiring a CR+LF command (*Carriage Return + Line Feed*), or select the options **+FF** (Form-Feed) for printers requiring a Form-Feed command added.



LOOK UP IN THE PRINTER MANUAL THE COMMANDS USED. IF THIS INFORMATION ARE NOT AVAILABLE, CARRY OUT A TEST PRINTOUT WITH EACH OF THE ABOVE OPTIONS TO INDIVIDUALIZE THE CORRECT OPERATING SETTING.

b) **REPORT** option

The following choice is proposed:

→	NR . COPIES	+
	PRINT LAST	↑
	EXIT	↓
		-

NR. COPIES allows to set the number of copies printed at the end of the sterilization process. Appears the following screen:

COPIES : XX	
+ / -	to set
↵	to enter
↑	to exit

Use the key + / - to set the desired number of copies (maximum 5). Confirm through the key ↵.

The item **PRINT LAST** allows to reprint the report of the last cycle performed (both correctly terminated or interrupted due to an alarm). The following choice is proposed:

→	NORMAL PRINT	+
	EXTENDED PRINT	↑
	EXIT	↓
		-

The option **NORMAL PRINT** enables a normal printout (that is a report with only the meaningful data of the cycle, and outputted at the end of a cycle correctly performed), whereas **EXTENDED** enables the complete printout (that is a report containing all the data, typical for a cycle interrupted by an alarm).



IF THE LAST CYCLE HAS BEEN CORRECTLY COMPLETED (OR INTERRUPTED THROUGH A MANUAL STOP) THE REPRINT IN NORMAL OR EXTENDED MODE CAN BE PERFORMED. IF THE LAST CYCLE HAS BEEN INTERRUPTED DUE TO AN ALARM OCCURRED (EXCLUDED A MANUAL STOP) ONLY THE **EXTENDED** MODE WILL BE AVAILABLE.

By entering the reprint command, the display will show the following message:

NOW PRINTING	
THE REPORT	
PLEASE WAIT . . .	

that will remain on until the end of the print.

Setting the tank Filling mode

(FILLING OPTIONS in ADVANCED menu)

Selecting **FILLING OPTIONS** the following choice is proposed:




millennium B_μ provides only the **MANUAL FILLING** option independently from the choice performed.

Setting the Draining mode

(DRAIN OPTIONS in ADVANCED menu)

Selecting **DRAIN OPTIONS** the following menu is shown:



millennium B_μ provides only the **INTERNAL DRAIN** option independently from the choice performed.

Acquisition of the environment pressure value

(option AMBIENT PRESSURE in menu SPECIAL)

At the first set-up and after each new installation of the sterilizer, the user must activate the acquisition of the environment pressure value. This operation is **necessary** for the correct operation of the internal auxiliary elements of the equipment.

By selecting the option **AMBIENT PRESSURE** the following screen will show:



CHECK THAT THE DOOR OF THE STERILIZER IS COMPLETELY OPEN. ON THE CONTRARY, IF YOU ATTEMPT TO PERFORM THE PRESSURE ACQUISITION WITH THE DOOR CLOSED, THE FOLLOWING MESSAGE WILL BE PRODUCED AND MAINTAINED UNTIL THE DOOR IS NOT OPEN:



Confirm by key ↓ to acquire the data. The following message will appear:



followed by an acoustic signal. The ambient pressure data is acquired.


Use key ↑ to exit this option without acquiring the data.

Adjusting the LCD contrast

(LCD CONTRAST option in SPECIAL menu)

Adjust the contrast of the LCD display to compensate different equipment positioning or environmental lighting conditions.

By selecting the option **LCD CONTRAST** the following choice is proposed:



Use key + or - to change the contrast accordingly.

Face the LCD display in normal working position and adjust the contrast for the better neatness of the writings.

Use the key ↑ to exit the option.

EXIT THE SETUP MODE


After having properly set any of the sterilizer options, proceed as follows to return in the normal operation mode:

- Move to main menu (see **Flowchart of the SETUP menu**).



TO RETURN IN MAIN MENU FROM WHATEVER PART OF THE PROGRAM YOU ARE, SIMPLY SELECT ONE OR MORE TIMES THE ITEM EXIT (CONFIRMING WITH KEY ↵)) DISPLAYED IN THE MENUS. ALTERNATIVELY PRESS THE KEY ⏏ (ESC) ONE OR MORE TIMES.

- Select the item **EXIT SETUP** and confirm with ↵. The display will show the message:



After a few seconds the equipment returns in normal operation mode and in **STAND-BY** status.

PREPARING THE MATERIAL TO BE STERILIZED

INTRODUCTION

The process sterilization can consider effective, reliable and repeatable on condition that the material is first suitably treated and subsequently tidy and correctly arranged into the sterilization chamber.

We notice that the organic residues or deposits of substances used in the medical practice are obviously receptacles of microorganisms and can hamper the contact of the steam with the instrument surfaces, inactivating, at least locally, the lethal process that the sterilization normally guarantees.

An incorrect arrangement of the load can lead to a difficult and sometimes impossible flowing and/or penetration of the steam on the material, with imaginable consequences. The drying process can strongly be also affected by this factor.

Therefore we suggest some **basic directions** concerning this aspects, leaving the user to deepen the problem in the most opportune way.

HANDLING THE MATERIAL BEFORE STERILIZING

Before everything we remember **some precautions** that are basic for the handling and moving the contaminated material:

- Wear rubber gloves of suitable thickness;
- Wash one's hands, already covered by the gloves, with a germicide detergent;
- Always use a tray to move the tools;
- Do never transport them directly picking up by hand;
- Protect the hands from the contact of possible sharpened or cutting parts in order to avoid the risk of dangerous infections;
- Immediately separate each item that has not to be sterilized or is not be able to support this process;
- Wash carefully one's hands covered with the gloves as the handling the not sterile material is over.

All the material and/or instruments to be sterilized should be then perfectly cleaned and deprived of whatever kind residues (organic and inorganic deposits, fragments of paper, buffer of cotton or gauze, calcareous residues etc.).



THE LACK OF THE CLEANING AND REMOVAL RESIDUE PROCEDURE, BESIDES TO CAUSE PROBLEMS DURING THE STERILIZATION PROCESS, CAN PROVOKE DAMAGES TO THE TOOLS AND/OR THE VERY SAME STERILIZER.

For an effective **cleaning**, proceed as follows:

1. Rinse the tools under a throw of running water, **immediately** after the use;
2. Separate the metallic instruments according to the material type (carbon steel, stainless steel, brass, aluminum, chrome, etc.) in order to avoid electrolyte oxidation;
3. Wash by using an ultrasonic equipment with a mixture of water and germicidal solution, taking care to follow the recommendations of the manufacturer.
4. For the best results use a neutral pH detergent specifically studied for the ultrasonic washing.



SOLUTIONS CONTAINING FENOLIS OR QUATERNARY AMMONIUM MIXTURES COULD CAUSE CORROSION ON THE TOOLS AND METALLIC PARTS OF THE ULTRASONIC EQUIPMENT.

5. After the washing, carefully rinse the tools and verify the complete elimination of the residues; if necessary repeat the washing cycle or operate manually.



IF POSSIBLE, USE DEJONIZED OR DISTILLED WATER FOR THE RINSING OPERATION IN ORDER TO PREVENT THE FORMATION OF CALCAREOUS STAINS. IF HIGH HARDNESS TAP WATER IS USED FOR THIS OPERATION, IT IS ADVISED TO ALWAYS DRY THE INSTRUMENTS.

For the **handpieces** (turbines, contrangles, etc.), integrate what above described with a treatment in special equipment for the internal cleansing and wiping (sometimes including the lubrication operation).



AT THE END OF THE STERILIZATION PROGRAM, REMEMBER TO LUBRICATE THE HANDPIECES INTERNAL MECHANISMS BY USING THE SPECIAL STERILE OIL. THIS PRECAUTION GUARANTEES THAT THE USEFUL LIFE OF THE TOOL DOESN'T RESULT IN ANY WAY REDUCED.





CONSULT THE INDICATIONS SUPPLIED BY THE MANUFACTURER OF THE INSTRUMENT/MATERIAL TO BE STERILIZED BEFORE TREATING IT INTO THE AUTOCLAVE, VERIFYING POSSIBLE INCOMPATIBILITIES. METICULOUSLY FOLLOW THE USE INSTRUCTIONS OF THE CLEANSING OR DISINFECTANT PRODUCTS AND THE OPERATING MANUAL OF THE WASHING AND/OR LUBRICATION AUTOMATIC EQUIPMENT .

As regards the textile material (or in general the porous material), as white uniforms, napkins, caps and other, provides for a careful wash and dry before treating it in autoclave.



DON'T USE DETERGENTS WITH HIGH CONTENT OF CHLORINE AND/OR PHOSPHATES. DON'T BLEACH WITH CHLORINE PRODUCTS. SUCH COMPONENTS COULD DAMAGE THE TRAY-HOLDER, TRAYS AND METALLIC TOOLS ARRANGED INTO THE STERILIZATION CHAMBER.

ARRANGEMENT OF THE LOAD



To get the better effectiveness of the sterilization process and to preserve the material in time, follow the indications below reported.

General notes for the arrangement of the load on the trays.


- Arrange the tools of different metal (stainless steel, moderate steel, aluminum, etc.) on different trays or however well separate between them.
- In case of **not** stainless steel tools, interpose a sterilization paper napkin or muslin cloth between tray and tool, avoiding direct contacts between the two different materials;
- Arrange however the objects sufficiently outdistanced, in order they maintain the layout for the whole sterilization cycle;
- Verify all the tools are sterilized in open position;
- Arrange the cutting tools, (scissors, lancets, etc.) so they cannot come in contact during the sterilization process; if necessary use a cotton cloth or gauze to isolate and protect them;
- Arrange the containers (glasses, cups, test-tubes, etc.) on one side or inverted position, avoiding possible water stagnation;
- **Don't overload the trays over the stated limit (see Appendix A)**. In a lot of situations, this maximum admitted value might be excessive; in this case always **use a bit of common sense**.
- **Don't stack** the trays one above the other or put them in direct contact with the walls of the sterilization chamber. Always use the **supplied tray-holder**.
- To introduce and extract the trays from the sterilization chamber, always use the supplied special tray removal tong.




TO DETECT THE PROCESS RESULT, SET A STERILIZATION INDICATOR TEST PER EACH TRAY: THIS PRECAUTION AVOIDS FROM PROCESSING ONCE AGAIN THE SAME LOAD OR, FOR THE WORSE, FROM USING A NOT STERILIZED MATERIAL. IF WRAPPED MATERIAL HAS TO BE TREATED, SET THE INDICATOR TEST INSIDE ONE OF THE WRAPS.

Indications for rubber and plastic pipes

- Always rinse with pyrogen-free water before the use; don't dry them;
- Arrange the pipes on the tray so that the extremities are not obstructed neither crushed.
- Don't provoke pleats neither winds, but leave stretched more linearly as possible.






Indications for wrapped material

- Wrap the tools one by one or, if more tools have to be set in the same wrap, verify that they are of the same metal;
- Seal the wrap with autoclave adhesive ribbon or by a thermal sealer.
- Don't use metallic clips, pins or other, as this jeopardizes the maintenance of the sterility;
- Arrange the envelopes in order to avoid the formation of air pockets that potentially could hamper the correct penetration and removal of the steam.
- Turn the envelopes in order to set the plastic part downward (tray side) and the paper part upward. Verify the correctness of this position, reversing it if necessary.
- If possible, arrange the envelopes on the edge by using a proper stand.
- Never pile up the envelopes.



ALWAYS WRAP THE TOOLS IN CASE OF PROLONGED STORE. SEE ALSO THE CHAPTER "PRESERVING THE STERILIZED MATERIAL" .




Use of the tray-holder, box-holder and instrument-holder box

The tray-holder and box-holder lock on a screw's head positioned at the end plate of the sterilization chamber; for unlock it lift slightly by hand the holder or the box-holder on the bottom of the chamber by using eventually a tool.



USE A PROPER GLOVE TO AVOID SCALDS.



The instrument-holder box is made of two parts locked by clips. Arrange the instruments as shown in the figure and the possible packets in the front space; then join by clips the two parts of the box and position it on the box-holder.

SELECTING THE STERILIZATION PROGRAM

INTRODUCTION


The program choice is fundamental for the success of the sterilization process.

Since every tool or material has conformation, consistence and different characteristics, it is very important to **individualize the more suitable program**, both to preserve its physical characteristics (avoiding or however limiting any alterations) and to guarantee the better sterilization effectiveness.



A GUIDE FOR SELECTING THE PROPER PROGRAM VERSUS THE LOAD IS SHOWN IN THE **APPENDIX B (PROGRAMMES)**.


ABOUT THE SELECTION



Switch on the equipment as described in the **Chapter "First start-up"**.



IF THE PASSWORD FEATURE HAS BEEN ENABLED WITH THE OPTION **ANY POWER ON OR ANY CYCLE START** (SEE CHAPTER PASSWORD **SETTING**), A SCREEN FOR THE PASSWORD ENTERING WILL BE SHOWN AT THE STERILIZER SWITCHING ON:



Enter the password previously set by using the keys + and -. Confirm by the key ↓.

The display doesn't propose any active pre-selection. The equipment is waiting for the program selection by the user.

Operate the **Program Selection** key, pressing it one or more times up to the desired program (preset program 1, 2, 3 or 4, properly signaled on the command panel).


LCD display shows on the higher lines the description of the selected cycle with the associated drying type, and on the lower lines the set-point values for temperature (**°C**), pressure (**bar**) and time (**mm:ss**) of the selected cycle.



ON PRESSING THE SELECTION KEY, THE FIRST PROPOSED STERILIZATION PROGRAM IS THE CYCLE LAST PERFORMED.

The display shows on the first two lines the program name and the drying mode set; on the next lines the temperature (**°C**), pressure (**bar**) and time (**mm:ss**) set-point values of the cycle.

For example, LCD will show as follows:



After a few seconds, the display switches to show the current temperature and pressure, and the current date and time.



To cancel the selection press shortly the ESC key (↑) on the command panel.



IF NO SELECTION IS ENTERED, THE EQUIPMENT CANNOT START ANYWAY. BY PRESSING THE START KEY WITHOUT HAVING SELECTED A PROGRAM, THE DISPLAY WILL SHOW THE FOLLOWING MESSAGE **ACCOMPANIED BY AN ACOUSTIC SIGNAL**:

S E L E C T A P R O G R A M
P L E A S E . . .



THE USE OF IMPROPER PROGRAM FOR THE TYPE OF MATERIAL TO BE STERILIZED (SEE APPENDIX B) **DOESN'T GUARANTEE** THE EFFECTIVENESS OF THE STERILIZATION PROCESS.

RUNNING THE STERILIZATION PROGRAM

GENERAL


The sterilization program runs along a cycle characterized by a lot of phases.

The number and the duration of the phases can differ between the PROGRAMMES depending on the air exhausting, sterilization process and drying modes.

The electronic control system monitors the running of the different phases, verifying at the same time if the parameters are correctly respected; if whatever type of anomaly is detected during the cycle, the program will immediately be interrupted and enter in alarm status identified by code and proper message showing the nature of the problem.

With this type of control, the perfect sterilization in every condition will be guaranteed whatever is the sterilization program selected.

STARTING THE PROGRAM



Now that the load is arranged into the sterilization chamber (with the precautions described on the **Chapter "Preparing the material to be sterilized"** and the desired program selected, **close** the door for an acoustic click.


The Door Status signaling  **is flashing** (door closed).

Push on the key **START**.
The door-locking mechanism is engaged.

The Door Status signaling  changes in **ON state** (door locked).




IF THE PASSWORD FEATURE HAS BEEN ENABLED WITH THE OPTION **ANY CYCLE START** (SEE CHAPTER PASSWORD **SETTING**), A SCREEN FOR THE PASSWORD ENTERING WILL BE SHOWN AT THE STERILIZER SWITCHING ON:



Enter the password previously set by using the keys **+** and **-**.
Confirm by the key **↵**.

After entered the **START** command and during all the sterilization cycle, the bottom lines of display will show the following values:



The counting of the time starts from the start command (first vacuum phase), excluding the preheating phase.

SEQUENCE OF THE PROCESS

Warm-up

It follows the description of the sterilization cycle, phase by phase. As example, we will use the most complete and meaningful cycle, i.e. the cycle relating the program **134°C POROUS** (preset 1 on the command panel), provided with fractionated pre-vacuum.

Entering the **START** command, the first phase performed is the preheating phase (**WARMUP**) in order to set the chamber temperature at the preset cycle beginning condition. The display is showing:



The process signaling  is off.


First vacuum pulse

Reached the stated heating conditions, starts the first vacuum phase (**1. VACUUM PULSE**) reducing the pressure into the chamber at the preset value.




First pressure pulse

Reached the preset vacuum value, the steam enters into the chamber and the pressure increases (**1. PRESSURE PULSE**) up to reach the preset value. The display shows:




Second vacuum pulse

At the end of the pressure phase the steam is discharged and starts the second vacuum phase of the sterilization chamber (**2. VACUUM PULSE**). The display shows.



Second pressure pulse


The second vacuum phase is followed by a new steam entry into the sterilization chamber, and consequently the pressure raises once again (**2. PRESSURE PULSE**).




The process signaling icon  is still switched off.

Third vacuum pulse


Follows a new discharge of the steam and the start of a third vacuum pulse (**3. VACUUM PULSE**). The display shows:

Third pressure pulse

The last vacuum phase is followed by the last steam entry, and the pressure increases once again (**3. PRESSURE PULSE**) up to the value preset for the sterilization process.


Thermodynamic equilibration


Reached the preset pressure and temperature values for the selected program, the program waits few seconds to allow the homogenization of the temperature inside the chamber and load (**EQUILIBRATION** time). The display is now showing:

Sterilization time

Stabilized the thermodynamic parameters, it begins the effective sterilization process of the material (**HOLDING TIME**). Thanks to the continuous monitor of the thermodynamic parameters and the full management of the hydraulic circuit, the pressure and temperature are maintained **constant** within the range stated by the program.

The countdown of the sterilization time starts .



Now the process signaling  **is blinking** to signal that the sterilization process is in progress.



OVER THE PROCESS PHASE, THE SIGNALING  CHANGES IN ON STEADY STATE TO INDICATE THE COMPLETION OF THE STERILIZATION.



IF, FOR WHATEVER REASON, THE STERILIZATION CYCLE IS INTERRUPTED **BEFORE** THE COMPLETION OF THE **PROCESS** PHASE, THIS SIGNALING REMAINS **BLINKING**. **IN THIS CASE THE MATERIAL INTO THE STERILIZATION CHAMBER CANNOT BE ANYWAY CONSIDERED STERILE AND MUST NOT ABSOLUTELY BE USED.**

Steam discharge


After the sterilization phase follows the steam exhausting from the sterilization chamber (**DEPRESSURIZATION**).




The process signaling  remains **ON**.

Drying


The natural depressurization is followed by a steam forced removal by means of the vacuum pump (**DRYING**) that generates in the chamber a additional depresseure to facilitate the steam exhausting. Now the LCD display one of the following screen depending on the drying mode set:



Standard drying




Intelligent drying



Extra drying
(+XX) means the
set time

Ventilation

The vacuum drying phase is followed by the ventilation phase (**VENTILATION**) during which, on maintaining the vacuum into the chamber, sterile and fresh air is entered to eliminate the condensation and cool the load.




Levelling

Ended the ventilation phase the chamber will be vented by letting the sterile air enter the chamber (**LEVELLING**) that will reach the atmospheric pressure value, and allow the door opening to recovery the load. LCD display changes in:




Cycle completion

As the pressure into the sterilization chamber reaches the pre-set safety limits, the door locking mechanism will be released.


Consequently the door status signaling  **starts to blink** and contemporary an acoustic warning signaling is generated.

Now the screen will show:



The process signaling  is still **ON**.

Open the door and recover the sterilized material by using the supplied tray removal tong.

Opening the door, the process signaling  **turns off**.



OPENING THE DOOR, A REPORT OF THE STERILIZATION CYCLE WILL BE PRINTOUT (IF A PRINTER IS INSTALLED AND SELECTED). CHECK THIS REPORT, SIGN AT THE PROPER LINE AND FILE IT IN A SURE PLACE. REFER TO THE REPORT EXAMPLES IN APPENDIX B, PROGRAMMES.


OPENING THE DOOR, THE EQUIPMENT GOES IN STAND-BY WITH THE MODE AS PREVIOUSLY SET. REPEAT THE PROCEDURES DESCRIBED IN THE CHAPTER "SELECTING THE STERILIZATION PROGRAM" IF A NEW STERILIZATION CYCLE HAS TO BE PERFORMED.

COMPLETED THE PROGRAM, AND UNTIL THE DOOR IS NOT OPEN, THE HEATING RESISTORS ARE DISABLED. ACCORDINGLY THE EQUIPMENT COOLS DOWN SLOWLY, WHATEVER IS THE CURRENT STAND-BY MODE.

Now the equipment is **ready** to perform a **new cycle**.






COMPLETED THE CYCLE, WHENEVER THE DOOR IS NOT OPEN, THE VACUUM PUMP IS PERIODICALLY ACTIVATED TO REMOVE ANY TRACE OF CONDENSATE FROM THE STERILIZATION CHAMBER. LCD IS SHOWING:



Press ↑ to stop the forced ventilation and open the door.


MANUAL INTERRUPTION OF THE PROGRAM



If necessary, the program can be manually interrupted by the operator in whatever moment by holding key **START/STOP** pressed for 3 seconds.


This command is detected by the equipment as an **alarm**, as the program is not correctly completed.

Accordingly the display is showing until the safety conditions are not reached:



with the error code **E999** on the display associated with a warning tone.

Reached the safety conditions, the machine activates a **special procedure** asking for a manual door unlocking by the user:



Push on the key **↑** to release the locking mechanism. The message changes in:



Finally, opening the door you are asked for **resetting** the equipment:



To reset the system **hold-down for 3 seconds and more the key PROGRAM SELECTION** until a confirmation tone signal.

Opening the door, the report of the sterilization cycle performed will be printout, including the error signaling (**E999**). Check this report, sign and file it in a sure place.

Refer to the report examples in **Appendix B, PROGRAMMES**.

After the RESET the equipment enters in STAND-BY state and now will be ready to perform a new cycle.




WHEN THE ALARM OCCURS DURING SOME CYCLE PHASES, AN AUTOMATIC HYDRAULIC CIRCUIT CLEANING PROCEDURE WILL BE ENABLED (SEE **APPENDIX E, ALARM SIGNALING**).



FOR THE COMPLETE DESCRIPTION OF THE ALARMS SEE **APPENDIX E (ALARM SIGNALING)**.



AFTER A MANUAL INTERRUPTION OF THE PROGRAM (MANUAL STOP) ALWAYS VERIFY THE STATE OF THE INDICATION  BEFORE USING THE TREATED MATERIAL OF THE STERILIZATION CHAMBER.

IF THE INDICATION  IS STEADY ON THE LOAD CAN BE CONSIDERED STERILE AND CONSEQUENTLY USED. WE RECOMMEND ITS IMMEDIATE USE.

ON THE CONTRARY, IF THE SIGNALING IS OFF, THE MATERIAL CONTAINED INTO THE STERILIZATION CHAMBER CANNOT BE CONSIDERED STERILE AND MUST NOT ABSOLUTELY BE USED.

RESULT OF THE PROGRAM

After the completion of the cycle we recommend to verify the result of the sterilization

If the cycle is normally over (indication **CYCLE COMPLETE**) without interruption caused by whatever alarms, the complete asepsis of the material will be guaranteed.

A further verification tool is represented by the printout of the sterilization parameters.

PRINTING THE DATA REPORT

(on external optional printer)

It is a good rule, to verify the positive result on the report issued at the end of the sterilization program.

The print of the meaningful data relating the thermodynamic parameters, pressure and temperature (bar and °C), and time (minutes) of the sterilization cycle, with special care on the sterilization phase, will be automatically performed at the end of every cycle simply by opening the door.

Check the values on the report and possible additional information to have further confirmation of the quality of the performed sterilization process.


Sign the report at the proper line and file it for a possible future use.

Copies of the report can be eventually used for identifying the load (or parts of it) with sterilization date/hour and details relative to the type of the performed cycle.

To choose the number of copies refer to the **Chapter "Setting up the equipment"**.



THE OPERATOR CAN SET TO GET AN EXTENDED REPORT OF THE STERILIZATION PROCESS DATA, INCLUDING ANY VALUES OF THE SENSORS INSTALLED ON THE EQUIPMENT. TO START THIS PRINTING OPTION HOLD DOWN THE KEY **ESC ↑** ON THE COMMAND PANEL WHILE OPENING THE DOOR.



For the details concerning the printing options, refer to the report examples in **Appendix B, Programmes**.

PRESERVING THE STERILIZED MATERIAL

GENERAL

The sterilized material has to adequately be handled and preserved to maintain own sterility in time until its use.

An improper maintenance **can** provoke its fast **recontamination**.

However this provokes an hazardous situation, as the alternative is to use the contaminated material (unconsciously in many cases) with risk for both operator and patient, or to perform a new sterilization process with inevitable waste of time and resources.

Therefore we retain useful to give some basic advice, letting the operator make further carefully investigations.

HANDLING

By assuming that the sterilizer is arranged in a clean place, without dust and damp, pay attention to the following indications when handling and moving the sterile material:

1. Remove the load from the sterilization chamber by wearing cleaned, or better, sterile gloves and white uniform. For greater precaution put a protective mask on the face;
2. Place the trays on a dry table, suitably cleaned and disinfected. Pay attention to outdistance or however separate the sterile material from the area of the contaminated material to be sterilized;
3. Touch the material and/or the tools as less is possible, paying very much attention do not lacerate or damage the wraps;
4. Let the tools cool before transporting (and storing). If necessary, use dry, clean and disinfected containers to move the material. The containers must be closed or, if open type, covered with clean cloths.

STORING

The sterile material, on waiting for the use, must be stored by taking opportune measures in order to slow as possible the contamination process:

1. Preserve the material and/or tools inside the protective wraps used for the sterilization. Don't wrap the tools after the sterilization as such practice, besides to be useless, is completely without meaning;
2. Store the material in a dry place, suitably clean and disinfected, away from areas where the infected material transits. If possible, prefer the closed areas provided with ultraviolet illumination;
3. Identify the sterile material by affixing the sterilization date (use a copy of the printed report or an adhesive labels);
4. In the first place use the material stored for longer time (by using a FIFO criteria, "first in first out"). This allows to have homogeneously stored material, avoiding too long storing periods with consequent risks.
5. Don't store for long time the material. Do not neglect that, even if the above indications have been followed, the material tends however to degrade, with a new contamination within a certain time.



CHECK ON THE SPECIFICATIONS SUPPLIED BY THE MANUFACTURER OF THE PACKING MATERIAL FOR THE MAXIMUM STORING TIME ADMITTED.

TEST PROGRAMMES

OVERVIEW

For the safety of both operator and patient, the functionality and effectiveness of an important process as the sterilization of medical devices must be periodically verified.

For this reason, **millennium B_μ** offers the possibility to perform in simple and automatic way two different test PROGRAMMES:

- Helix/BD Test
- Vacuum Test

Helix/BD Test allows to perform a cycle at 134°C characterized by a special time of the sterilization phase (3.5 min.); the cycle is provided with fractionated vacuum like the POROUS and HOLLOW PROGRAMMES.


Through a proper device it is possible to value the correct penetration of the steam inside the porous loads (see next paragraph).

Such cycle is also proper to measure the penetration of the steam inside the porous loads (standard **Bowie & Dicks** test pack).


Vacuum Test allows to verify the perfect watertight of the hydraulic circuit of the sterilizer.

By measuring the vacuum offset within a stated time and comparing it with preset limit values, it is possible to determine the quality of the sterilization chamber and connecting pipelines watertight.

HELIX/BD TEST



To select the Helix/BD Test program press one or two times the **TEST SELECTION** key until the display is showing:



The test device (according with the requirements of EN 867-5) is a PTFE pipe, long 1.5 m and of internal diameter of 2 mm, with a small watertight screw capsule mounted on one end containing an opportune chemical indicator.

The other end of the pipe is free to allow the penetration of the steam and value the effectiveness.

To perform the test (as stated in EN 867-5) insert the chemical indicator, consisting of a paper strip with a special reacting inch, inside the capsule of the device (that must be always perfectly dried). Close the capsule in such a way no leakage through the gasket is possible.




THE DEVICE AND THE CHEMICAL INDICATORS FOR THE HELIX/BD TEST ARE NOT SUPPLIED WITH THE EQUIPMENT. FOR MORE INFORMATION CONTACT "ASSISTENZA CLIENTI M.O.COM." (SEE APPENDIX Z)

Set the device on the middle tray of the equipment, approximately in the center. **Don't introduce** other material inside the chamber. Close the door and start the program by key **START**.



IF THE PASSWORD FEATURE IS ENABLED WITH THE OPTION **EACH CYCLE START** (SEE CHAPTER PASSWORD **SETTING**), A SCREEN FOR THE PASSWORD ENTERING WILL BE SHOWN ON SELECTING THE PROGRAM:



Enter the password previously set by using the keys + and -. Confirm by the key ↵.



ALL THE INDICATIONS DISPLAYED DURING THE CYCLE ARE THE SAME AS DESCRIBED IN CHAPTER "RUNNING A STERILIZATION PROGRAM".

At the end of the program, recover the test device from the sterilization chamber, open the capsule and remove the indicator from its place.

If the steam is correctly penetrated, the ink color on the whole length of the strip will be completely turned compared to the start conditions; on the contrary (insufficient penetration) the change will be only partial or missing.



NORMALLY THE TURNING OCCURS FROM A CLEAR COLOR (BEIGE, YELLOW, ETC.) TOWARD A DARK COLOR (BLUE, VIOLET OR BLACK). IN EVERY CASE, METICULOUSLY FOLLOW THE INSTRUCTIONS SUPPLIED BY THE MANUFACTURER ABOUT THE WAY FOR USING AND INTERPRETING THE INDICATOR DEVICE.

The test time takes about **15 minutes**.

Opening the door at the end of the test, and if a printer is installed, a report of the main data of the performed test will be released.

Affix the chemical indicator in the available area, sign the report and file it.

For more details about the printing reports, refer to the examples in **Appendix B, PROGRAMMES**.

VACUUM TEST



To select the VACUUM TEST press one or more times the SELECTION TEST key until the display is showing:




The Vacuum Test program must be carried out with the sterilization chamber empty, only trays and tray-holder should be loaded.



CARRY OUT THE VACUUM TEST AS FIRST CYCLE AT THE SWITCHING ON OF THE EQUIPMENT.

IN ORDER TO AVOID THAT THE HEATING EFFECT INFLUENCES THE VACUUM CHANGE VALUE READ DURING THE VACUUM TEST, THE SYSTEM IS PROGRAMMED TO DISABLE THE TEST IF THE CHAMBER TEMPERATURE SENSOR IS READING A VALUE HIGHER THAN 50° C.

If you are trying to launch the test program with chamber temperature higher than the above preset value, LCD will show the following:



After a few seconds, the unit will automatically return in STAND-BY mode, and will be ready for a new command.




TO ALLOW THE CHAMBER TEMPERATURE TO QUICKLY FALL, LEFT OPEN THE DOOR UNTIL A LOWER TEMPERATURE IS REACHED.

Close the door and launch the program through the **START** key.




IF THE PASSWORD FEATURE IS ENABLED WITH THE OPTION **EACH CYCLE START** (SEE CHAPTER PASSWORD **SETTING**), A SCREEN FOR THE PASSWORD ENTERING WILL BE SHOWN ON SELECTING THE PROGRAM:



Enter the password previously set by using the keys + and -. Confirm by the key ↓.


The vacuum phase immediately starts and the display will show:



The display is showing the pressure value (bar) and the total time from the start of the cycle.


As soon as the preset pressure value (**-0.80 bar**) is reached the vacuum pump is arrested and the pressure stabilization phase (**WAITING PERIOD**) enabled for 5 minutes.

The display is showing the following message:



During this phase an offset of only 10% of the vacuum maximum value is admitted without bearing a test failure. The time is counted down up to the completion of the phase.

Expired this time the effective pressure detection phase (**LEAKAGE PERIOD**) starts with a duration of 10 minutes. The display will show:




During this phase an offset of only ± 0.02 bar from the start value is admitted. Higher offsets lead to the failure of the test.

The time is counted down up to the completion of the phase.

Expired this time the pressure is returned to the atmospheric value and the display will show:



Over the vacuum test, the display is showing:



A signaling tone points out the end of the program.



IF THE PROGRAM IS INTERRUPTED FOR A PRESSURE OFFSETS EXCEEDING THE PRE-SET LIMITS, AN ALARM MESSAGE WILL BE PRODUCED. FOR THE DESCRIPTION OF THE ALARMS SEE APPENDIX E , **SOLUTION OF THE PROBLEMS** .

The test takes about **16 minutes**.

If a printer is installed, on opening the door at the end of the program a report of the performed test will be released, including the main data.

For more details about the printing reports, refer to the examples in **Appendix B, PROGRAMMES**.

APPENDIX A – TECHNICAL CHARACTERISTICS

OVERVIEW TABLE

Equipment	Steam sterilizer	
Classification (as per 93/42/CEE)	Ila	
Model	millennium[®] B_μ	
Manufacturer	M.O.COM. S.r.l. Via delle Azalee, 1 20090 BUCCINASCO (MI) - ITALY	
Power supply	220V - 240 V~	
Frequency	50/60 Hz	
Mains fuses (6,3 x 32 mm)	F 16A 250V	
Power supply board fuses (5 x 20 mm)	F1: T 5A 250V (secondary winding of the trafo – 24V) F2: T 2A 250V (5V digital circuit / 12V analogue circuit) F3: T 2A 250V (primary winding of the trafo – 230V) F4: F 200mA 250V (door-locking protection) F5: F 1.25A 250V (door-locking motor reducer overload) F1 PTR: T 4 A (printer protection)	
External dimensions (LxDxH)	380 x 585 x 340 mm (rear connections excluded)	
Nominal power	3000 W (13 A)	
Insulation class	Class I	
Installation category	Cat. II	
Utilization	Internal use	
Acoustic noise level	60 db(A) max	
Environmental operating conditions	Temperature: +15°C ÷ +40°C Relative Humidity: max 80% , non condensing Height: max 3000 m (a.s.l.)	
Net weight	approx. 38 kg (empty) approx. 40 kg (empty, with trays and tray-holder) approx. 44 kg (empty, with trays + tray-holder + water at Max level)	
Sterilization chamber dimensions (L _{max} x H _{max} x D)	215 x 76 x 305 mm	
Total volume of the sterilization chamber	approx. 5.5 l	(0.0055 m ³)
Useful volume of the chamber (with tray-holder placed)	approx. 3 l	(0.003 m ³)
Capacity of the distilled water reservoir (feeding)	approx. 3.6 l approx. 0.8 l	(water at level MAX) (water at level MIN)
Sterilization programmes	Available: Preset:	11 (see Appendix B) 4 (directly selected by the operator)
Test programmes	Helix/BD Test Vacuum Test	
Pre-heating time (from cold)	approx. 10 min	
Printer parallel interface	DB-25 pin connector	(female)
PC serial interface	DB-9 pin connector	(female)
Bacteriologic filter (PTFE)	Porosity: Connection:	0,2 μm 1/8" NPT male connector

SAFETY DEVICES

millennium B_μ is provided with the followings safety devices; a brief description of their function is given on the following:

- **Mains fuses** (see the data on the overview table)
Protection of the whole equipment against possible failures of the heating resistors.
Action: interruption of the electric power supply.
- **Protection fuses on electronic circuits** (see the data on the overview table)
Protection against possible failures of the mains transformer primary circuit and low voltage loads.
Action: interruption of one or more low voltage circuits.
- **Thermal cutouts on the mains transformer windings**
Protection against possible vacuum pump and mains transformer primary winding overheating.
Action: temporary interruption (up to the cooling) of the winding.
- **Safety valve**
Protection against possible sterilization chamber over-pressure.
Action: release of the steam and restoration of the safety pressure.
- **Manually resettable thermostat on steam generator**
Protection for possible steam generator overheating.
Action: interruption of the steam generator power supply.
- **Manually resettable thermostat on chamber heating resistors**
Protection for possible overheating of the chamber heating resistors.
Action: interruption of the power supply of the chamber resistors.
- **Safety micro-switch for the door status**
Comparison for the correct closing position of the sterilization chamber door.
Action: signaling of wrong position of the door.
- **Door locking motorized mechanism with electromechanical protection (pressure switch)**
Protection against accidental opening of the door (also in case of black-out).
Action: impediment of the accidental opening of the door during the program.
- **Safety micro-switch on the door locking mechanism**
Comparison for the correct closing position of the locking system.
Action: signaling of the unsuccessful or incorrect operation of the door locking mechanism.
- **Self-leveling hydraulic system**
Hydraulic system for the natural pressure leveling in case of manual cycle interruption, alarm or black-out.
Action: automatic restoration of the atmospheric pressure inside the sterilization chamber.
- **On-board system for the sterilization process evaluation**
Continuous monitor of the sterilization process parameters, fully microprocessor managed.
Action: immediate interruption of the program (in case of anomaly) and generation of alarms.
- **Monitoring of the sterilizer operation**
Real time monitoring of all the meaningful equipment parameters.
Action: generation of alarm messages (in case of anomaly) with possible interruption of the cycle.

CHARACTERISTICS OF THE FEEDING WATER

DESCRIPTION	VALUES IN THE DISTILLED WATER	VALUES IN THE CONDENSED
DRY RESIDUE	< 10 mg/l	< 1 mg/l
SILICON MONOXIDE SiO ₂	< 1 mg/l	< 0,1 mg/l
IRON	< 0,2 mg/l	< 0,1 mg/l
CADMIUM	< 0,005 mg/l	< 0,005 mg/l
LEAD	< 0,05 mg/l	< 0,05 mg/l
HEAVY METAL RESIDUES (except iron, cadmium and lead)	< 0,1 mg/l	< 0,1 mg/l
CHORIDES	< 2 mg/l	< 0,1 mg/l
PHOSPHATES	< 0,5 mg/l	< 0,1 mg/l
CONDUCTIVITY AT 20°C	< 15 µs/cm	< 3 µs/cm
VALUE of pH	5 - 7	5 - 7
ASPECT	colorless, transparent, without sediments	colorless, transparent, without sediments
HARDNESS	< 0,02 mmol/l	< 0,02 mmol/l



ON BUYING THE DISTILLED WATER, ALWAYS CHECK THE QUALITY AND THE CHARACTERISTICS DECLARED BY THE MANUFACTURER ARE COMPATIBLE WITH THE ONES REPORTED IN THE ABOVE TABLE.



THE USE OF STEAM GENERATED WATER WITH CONTAMINANTS EXCEEDING THE VALUES INDICATED ON THE ABOVE TABLE CAN NOTABLY SHORTEN THE LIFE OF THE STERILIZER. BESIDES THIS CAN PRODUCE AN INCREASE OF THE OXIDATION ON THE MOST SENSITIVE MATERIALS AND AN INCREASE OF THE CALCAREOUS RESIDUES ON GENERATOR, CHAMBER, INTERNAL SUPPORTS, TRAYS AND TOOLS.

APPENDIX B – PROGRAMMES

INTRODUCTION

The steam sterilization is suggested for nearly all the materials and tools on condition that these are able to support without damages a minimum sterilization temperature of **121°C** (on the contrary other sterilization systems at lower temperature should be used).

The material normally autoclavable into a steam sterilizer is as follows:

- Stainless steel surgical instruments;
- Stainless steel generic tools;
- Carbon steel generic tools;
- Dynamic instruments, motorized by air (turbines) or gears (contrangles, ablation tools, etc.);
- Glass articles;
- Mineral basis articles;
- Heat-proof plastic articles;
- Heat-proof rubber articles;
- Heat-proof textile articles;
- Treatment material (gauze, tampons, etc.);
- Other autoclavable generic material.



DEPENDING ON THE CONFORMATION (SOLID, HOLLOW OR POROUS), PACKING (PAPER/PLASTIC WRAPPING, STERILIZATION PAPER, CONTAINER, MUSLIN ETC.) AND HEAT-PROOF CHARACTERISTIC OF THE MATERIAL, THE CORRECT STERILIZATION PROGRAM MUST BE SELECTED BY REFERRING TO THE TABLE OF THE NEXT PAGE.



THE RAPID PROCESS PROVIDED BY THIS EQUIPMENT MAKES INADVISABLE THE USE OF CLOSED CONTAINERS FOR THE STERILIZATION. PLEASE CONTACT "ASSISTENZA CLIENTI MOCOM" (SEE APPENDIX Z) FOR MORE DETAILS.



THE DEVICE MUST NOT BE USED FOR THE STERILIZATION OF FLUIDS, LIQUIDS OR PHARMACEUTICAL PRODUCTS.



OVERVIEW OF THE AVAILABLE PROGRAMMES



PROGRAM DESCRIPTION	NOMINAL VALUES				BASIC PARAMETERS OF THE PROGRAM					STERILIZABLE MATERIAL				NOTES
	Temperature (°C)	Pressure (bar)	Holding time (min)	Cycle type (prEN 13060: 2004)	Pre-vacuum (F= fractionated; S= single)	Standard drying (L= long; S= short)	Total time of the cycle (medium load + max load)	Average water consumption of the cycle	Average power consumption (kWh / cycle)	TYPE	MAX LOAD TOTAL (kg)	MAX LOAD PER TRAY (kg)	MAX LOAD PER ITEM (kg)	
134°C POROUS	134	2,10	4	B	F	L	17+19	200	0,6	Unwrapped porous material	0.30	0.30	0.30	For porous wrapped or unwrapped material and (single or dual) wrapped instrument we suggest to use 1-tray only configuration
										Single-wrapped porous material	0.30	0.30	0.30	
										Dual-wrapped porous material	0.30	0.30	0.30	
										Single-wrapped solid and hollow material	1.00	0.50	0.25	
										Dual-wrapped solid and hollow material	0.50	0.25	0.25	
134°C PRION	134	2,10	18 >	B	F	L	31+33	210	0,7	Unwrapped porous material	0.30	0.30	0.30	
										Single-wrapped porous material	0.30	0.30	0.30	
										Dual-wrapped porous material	0.30	0.30	0.30	
										Single-wrapped hollow material	1.00	0.50	0.25	
										Dual-wrapped solid and hollow material	0.50	0.25	0.25	
121°C POROUS	121	1,10	20	B	F	L	32+34	210	0,6	Unwrapped porous material	0.30	0.30	0.30	
										Single-wrapped porous material	0.30	0.30	0.30	
										Dual-wrapped porous material	0.30	0.30	0.30	
										Single-wrapped hollow material	1.00	0.50	0.25	
										Dual-wrapped solid and hollow material	0.50	0.25	0.25	
134°C HOLLOW	134	2,10	4	S	F	S	15+17	200	0,6	Unwrapped hollow material	1.80	0.90	0.25	
121°C HOLLOW	121	1,10	20	S	F	S	30+32	210	0,6	Unwrapped hollow material	1.80	0.90	0.25	
134°C WRAPPED	134	2,10	4	S	S	L	14+16	125	0,5	Single-wrapped solid material	1.00	0.50	0.25	We suggest to use 1-tray only configuration
121°C WRAPPED	121	1,10	20	S	S	L	29+31	130	0,5	Single-wrapped solid material	1.00	0.50	0.25	
134°C SOLID	134	2,10	4	N	S	S	11+13	125	0,4	Unwrapped solid material	1.80	0.90	0.25	
121°C SOLID	121	1,10	20	N	S	S	26+28	130	0,4	Unwrapped solid material	1.80	0.90	0.25	
134°C EMERGENCY	134	2,10	3	N	S	Fast	9	125	0,4	Unwrapped solid material	1.80	0.90	0.25	
XXX°C USER (see note)	134 or 121	2,10 or 1,10	4 > or 20 >	n.a.	F/S	L/S	n.a.	n.a.	n.a.	Material and instruments according with the user setup	n.a.	n.a.	n.a.	Variable parameters depending on the user setup
HELIX/BD TEST	134	2,10	3,5	-	F	S	15	-	-	Test pack only (without any load)	-	-	-	
VACUUM TEST	-	-0,80	-	-	-	-	16	-	-	Chamber empty	-	-	-	







NO PASSWORD IS REQUIRED TO ACCESS THE USER PROGRAMME. ANY POSSIBLE PARAMETER COMBINATIONS OF THE CYCLE CUSTOMIZATION DO NOT INVOLVE AZARDS TO USER OR EQUIPMENT.



DIAGRAMS OF THE STERILIZATION PROGRAMMES













DIAGRAMS OF THE TEST PROGRAMMES



EXAMPLES OF PRINTING REPORTS

Normal program report

Model MILLENNIUM Bp
S/N 02 MM 0001
Ver. SW Exxxx/BPyyyyyy
Counter 0007/0015
Selection 134c SOLID
Temperature 134 °C
Pressure 2.10 bar
Process time 4 min
Stand-by LOW
Pre-vacuum SINGLE
Drying FAST

CYCLE START 19/01/04
12:14

Time		C	bar
00:01	CS	079.4	+0.00
01:45	1PV	093.7	-0.80
03:58	ET	135.6	+2.15
04:02	SS	135.9	+2.17
04:52		135.6	+2.14
05:40		135.5	+2.14
06:32		135.4	+2.14
07:30	SE	135.5	+2.15
07:55	DS	104.1	+0.00
08:50	SPD	047.5	-0.90
11:55	DE	047.6	-0.84
12:50	CE	084.6	-0.04

04:25 MAX 136.0
06:30 MIN 135.4

Drying Pulses 01
CYCLE END 19/01/04
12:23

STERILIZATION: POSITIVE

OPERATOR

Model MILLENNIUM Bp
S/N 02 MM 0001
Ver. SW Exxxx/BPyyyyyy
Counter 0007/0015
Selection 134c POROUS
Temperature 134 °C
Pressure 2.10 bar
Process time 4 min
Stand-by HIGH
Pre-vacuum FRACTIONATED
Drying STANDARD

CYCLE START 19/01/04
09:52

Time		C	bar
00:01	CS	075.1	-0.00
01:22	1PV	047.5	-0.80
03:25	1PP	120.5	+1.00
04:30	2PV	061.1	-0.80
06:05	2PP	120.4	+0.98
07:40	3PV	061.1	-0.80
10:12	ET	135.5	+2.15
10:25	SS	135.9	+2.17
11:09		135.4	+2.14
11:48		135.5	+2.15
12:04		135.4	+2.14
12:55	SE	135.5	+2.15
13:40	DS	104.4	+0.00
14:35	SPD	048.4	-0.90
18:55	EPD	094.9	-0.86
20:35	DE	112.6	-0.47
21:02	CE	115.8	-0.04

10:35 MAX 135.9
12:11 MIN 135.4

Drying Pulses 05
CYCLE END 19/11/02
10:09

STERILIZATION: POSITIVE

OPERATOR

Extended program report (required by operator)

Model MILLENNIUM Bp
S/N 02 MM 0001
Ver. SW Exxxx/BPyyyyyy
Counter 0007/0015
Selection 134c POROUS
Temperature 134 °C
Pressure 2.10 Bar
Process time 4 min
Stand-by HIGH
Pre-vacuum FRACTIONATED
Drying STANDARD

CYCLE START 19/01/04
09:52

Time		T1	P	T2	T3	T4
00:01	CS	075.1	-0.00	130.9	115.2	093.4
00:08		074.9	-0.28	133.3	114.2	094.0
00:15		074.4	-0.46	146.3	113.2	094.5
00:22		074.3	-0.57	152.6	112.2	095.0
00:28	..	074.3	-0.59	154.2	111.9	095.2
00:36	..	078.9	-0.62	152.2	110.4	095.6
00:45	..	074.9	-0.73	146.6	109.6	095.7
01:03	..	047.8	-0.78	149.3	107.7	095.7
01:25	..	047.8	-0.80	155.3	105.8	095.4

01:35 .. 076.5 -0.57 149.9 105.2 095.1
01:45 .. 081.1 -0.49 142.1 104.6 094.6

05:45 ... 068.4 -0.76 151.8 104.7 102.3
05:52 ... 061.1 -0.80 153.6 104.5 101.7

06:02 ... 097.4 +0.01 154.7 104.0 100.8
06:12 ... 104.6 +0.24 148.9 103.7 101.0

10:33 ... 135.5 +2.15 143.3 111.7 131.7

10:45 ... 135.9 +2.17 148.5 113.5 132.6

10:55 ... 135.3 +2.16 153.6 115.9 133.0

13:35 ... 135.5 +2.15 157.4 126.5 132.5

13:45 ... 134.4 +1.07 157.0 126.8 131.2

13:58 ... 108.3 +0.25 156.4 126.8 119.9

14:03 .. 104.4 +0.00 156.1 126.6 116.2

14:13 ... 094.2 -0.50 155.1 125.9 112.4

14:23 ... 069.2 -0.73 153.7 124.5 112.9

14:37 ... 059.2 -0.81 152.3 123.4 113.5

14:50 ... 053.8 -0.87 151.2 122.9 113.6

15:02 ... 048.4 -0.90 150.9 122.7 113.5

15:10 ... 047.1 -0.80 151.0 122.5 113.5

16:31 ... 042.3 -0.89 153.3 122.0 112.2

19:01 ... 094.9 -0.90 153.3 121.7 112.3

19:12 ... 101.4 -0.67 154.0 121.7 112.3

19:25 ... 105.4 -0.57 153.7 121.5 112.3

20:28 ... 112.6 -0.47 149.6 119.1 111.2

20:39 ... 115.2 -0.10 143.0 118.4 110.7

20:55 CE 115.8 -0.04 147.4 110.1 110.7

11:20 MAX 135.9

13:11 MIN 135.4

Drying pulses 05

CYCLE END 19/01/04
10:09

STERILIZATION: POSITIVE

OPERATOR

EXTENDED REPORT
REQUESTED BY THE OPERATOR

Report following a Manual Stop

Model MILLENNIUM Bp
S/N 02 MM 0001
Ver. SW Exxxx/BPyyyyyy
Counter 0007/0015
Selection 134c POROUS
Temperature 134 °C
Pressure 2.10 bar
Process time 4 min
Stand-by HIGH
Pre-vacuum FRACTIONATED
Drying STANDARD

CYCLE START 19/01/04
11:13

Time		C	bar
00:01	CS	077.6	+0.01
01:15	1PV	088.7	-0.80
03:16	1PP	120.6	+1.00
04:02	2PV	062.9	-0.80
05:10	2PP	135.6	+1.00
05:55	3PV	135.5	-0.80
08:02	ET	135.4	+2.15
08:17	SS	135.5	+2.17
09:02		135.5	+2.14
09:42		104.1	+2.15
10:39		047.5	+2.15

STERILIZATION: NEGATIVE

OPERATOR

ALARM CODE: E999
DESCRIPTION MANUAL STOP

Report following a Black-out

Model MILLENNIUM Bp
S/N 02 MM 0001
Ver. SW Exxxx/BPyyyyyy
Counter 0006/0012
Selection 134c CUSTOM
Temperature 134 °C
Pressure 2.10 bar
Process time 07 min
Stand-by HIGH
Pre-vacuum FRACTIONATED
Drying FAST

CYCLE START 19/01/04
15:31

BLACK OUT 19/01/04
15:39

STERILIZATION NEGATIVE

OPERATOR

ALARM CODE: E000
DESCRIPTION BLACK-OUT

Report following an Alarm

Model MILLENNIUM B μ
 S/N 02 MM 0001
 Ver. SW Exxxx/BPyyyyyy
 Counter 0007-0015
 Selection 134c POROUS
 Temperature 134 °C
 Pressure 2.10 Bar
 Process time 4 min
 Stand-by HIGH
 Pre-vacuum FRACTIONATED
 Drying STANDARD

CYCLE START 19/01/04 11:30

Time	T1	P	T2	T3	T4
00:01 CS	075.1	-0.00	130.9	115.2	093.4
00:08	074.9	-0.28	133.3	114.2	094.0
00:16	074.4	-0.46	146.3	113.2	094.5
00:23	074.3	-0.57	152.6	112.2	095.0
00:27 ..	074.3	-0.59	154.2	111.9	095.2
00:39 ..	078.9	-0.62	152.2	110.4	095.6
00:46 ..	074.9	-0.73	146.6	109.6	095.7
01:01 ..	047.8	-0.78	149.3	107.7	095.7
01:23 ..	047.8	-0.80	155.3	105.8	095.4
01:30 ..	076.5	-0.57	149.9	105.2	095.1
01:40 ..	081.1	-0.49	142.1	104.6	094.6
05:46 ...	068.4	-0.76	151.8	104.7	102.3
05:53 ...	061.1	-0.80	153.6	104.5	101.7
06:01 ...	097.4	+0.01	154.7	104.0	100.8
06:11 ...	104.6	+0.24	148.9	103.7	101.0
10:38 ...	135.5	+2.15	143.3	111.7	131.7
10:50 ...	135.9	+2.17	148.5	113.5	132.6
10:58 ...	135.3	+2.16	153.6	115.9	133.0
13:34 ...	135.5	+2.15	157.4	126.5	132.5
13:48 ...	134.4	+1.07	157.0	126.8	131.2
13:59 ...	108.3	+0.25	156.4	126.8	119.9
14:04 DS	104.4	+0.00	156.1	126.6	116.2

STERILISATION NEGATIVE

ALARM CODE: A112
 DESCRIPTION PTC SHORTCIRCUIT

CAUTION !
 PLEASE REFER TO USER MANUAL

HELIX/BD TEST program report

Model MILLENNIUM B μ
 S/N 02 MM 0001
 Ver. SW Exxxx/BPyyyyyy
 Counter 0011/0019
 Selection HELIX TEST
 Temperature 134 °C
 Pressure 2.10 bar
 Process time 3.5 min

CYCLE START 19/01/04 16:38

Time	C	bar
00:01 CS	076.4	+0.00
01:35 1PV	089.3	-0.89
03:20 1PP	120.4	+0.99
04:12 2PV	062.5	-0.78
05:02 2PP	120.2	+0.97
05:50 3PV	061.1	-0.79
07:45 ..	135.6	+2.15
07:55 ..	136.0	+2.17
08:40 ..	135.6	+2.14
09:25 ..	135.6	+2.15
10:12 ..	135.5	+2.14
10:34 ..	135.4	+2.14
10:50 ..	111.5	+0.00
11:40 ...	047.8	-0.89
12:55 ...	059.5	-0.86
13:45 ..	075.4	-0.50
14:34 CE	078.7	-0.04
08:05 MAX	136.0	
09:44 MIN	135.4	

Drying pulses 01
 CYCLE END 19/01/04 16:49

HELIX TEST COMPLETE
 Please attach the indicator hereunder

OPERATOR

VACUUM TEST program report

Model MILLENNIUM B μ
 S/N 02 MM 0001
 Ver. SW Exxxx/BPyyyyyy
 Counter 0011/0019
 Selection VACUUM TEST

CYCLE START 19/01/04 11:37

Time	C	bar
00:00 CS	035.0	+0.00
01:09 E1F	037.4	-0.80
6:09 E2F	038.4	-0.79
16:09 E3F	042.0	-0.79
17:04 CE	045.5	-0.01
CYCLE END	19/11/02 11:41	
VACUUM TEST:	POSITIVE	
OPERATOR		

APPENDIX C - MAINTENANCE

In order to guarantee a sure and efficient operation for the whole life of the equipment, a regular maintenance by the user as well as the correct use are necessary.

GENERAL

For a better quality of the maintenance, it is necessary to integrate the ordinary controls with periodic check-up by the Service (refer to appendix Z).

In addition, a **periodic validation of the sterilizer**, i.e. the check of the thermodynamic process parameters comparing them to the reference values measured by tools suitably adjusted, is fundamental.

Refer to the paragraph **Validation of the sterilizer**, in the following of this appendix.

The ordinary maintenance, as described below, consists of easy manual operations and preventive interventions by the use of simple tools.



IN CASE OF REPLACEMENT OF COMPONENTS OR EQUIPMENT PARTS REQUIRE FOR AND/OR USE ONLY ORIGINAL SPARE PARTS.

ORDINARY SCHEDULED MAINTENANCE


The following table reports an overview of the maintenance actions and frequency to be performed on the sterilizer in order to maintain it always efficient.

We suggest to shorten the maintenance interval in case of **heavy use**:

DAILY	Cleaning of the gasket and chamber cover. Cleaning of the external surfaces
WEEKLY	Cleaning of the sterilization chamber and accessories Desinfection of the external surfaces
MONTHLY	Lubrication of the door locking mechanism Maintenance of the safety valve Cleaning (or replacement) of the draining filter Sterilization of the bacteriologic filter
EVERY 3/6 MONTHS (depending on the use frequency)	Replacement of the bacteriologic filter
YEARLY	Validation of the sterilizer (<i>see proper paragraph</i>)

Always refer to the **following general directions**:

- **Do not wash** the sterilizer with direct pressurized or rain throw water. Possible water penetrations on electric and electronic components could jeopardize, also irreparably, the operation of the equipment or internal parts;
- **Do not use** metal cleaning abrasive cloths, metallic brushes (or other abrasive materials) or both solid and liquid products for cleaning the equipment or the sterilization chamber;
- **Do not use** chemical products or disinfectant substances for cleaning the sterilization chamber. These products could provoke possible damages, sometimes irremediable, to the sterilization chamber;
- **Do not let** the calcareous or dirty residues accumulate on the sterilization chamber, door and gasket, but provide for their periodic removal. In the course of the time these residues could cause damages to the elements, and jeopardize the operation of the components installed on the hydraulic circuit.



THE PRESENCE OF WHITE SPOTS AT THE BASE OF THE INSIDE WALL OF THE CHAMBER MEANS THE USE OF POOR QUALITY DEMINERALIZED WATER.



BEFORE CARRYING OUT THE ORDINARY MAINTENANCE OPERATIONS, CHECK THAT THE MAINS SWITCH OF THE STERILIZER IS IN OFF POSITION (EQUIPMENT SWITCHED OFF) . IF THE DISCONNECTION OF THE ELECTRIC POWER SUPPLY OF THE EQUIPMENT IS NOT POSSIBLE AND THE EXTERNAL MAINS BREAKER IS FAR OR HOWEVER NOT VISIBLE TO THE PERSON IN CHARGE OF THE MAINTENANCE, THE NOTICE WORK IN PROGRESS MUST BE POSITIONED ON THE BREAKER AFTER HAVING SWITCHED IT IN OFF POSITION.

MAINTENANCE ACTIVITIES

Cleaning the gasket and inner part of the door

With reference to the previous table, the operations concerning the different maintenance actions are here synthetically described.

In order to eliminate the calcareous residues, clean the inner part of the door and gasket with a cotton cloth soaked with bland water and vinegar solution (or analogous products, verifying in advance the content reported on the label).
Dry the surfaces and remove every possible residue before using the equipment.

Cleaning the external surfaces

Clean all the external parts by using a clean cotton cloth dampened with water, eventually added with some neutral detergent.
Dry the surfaces and remove every possible residue before using the equipment.

Cleaning the sterilization chamber and accessories

Clean sterilization chamber, tray-holder and trays (and all the internal surfaces) with a clean cotton cloth soaked with water, eventually added with a bland neutral detergent.
Rinse carefully with distilled water, paying attention do not let any residue remain on the chamber or accessories.



DO NOT USE POINTED OR SHARP TOOLS TO REMOVE THE CALCAREOUS SPOTS FROM THE STERILIZATION CHAMBER.

IN CASE OF EVIDENT RESIDUES, IMMEDIATELY VERIFY THE QUALITY OF THE DISTILLED WATER YOU ARE USING (SEE [APPENDIX A](#), TECHNICAL CHARACTERISTICS).

Disinfecting the external surfaces

For the desultory disinfection of the external surfaces both methylated spirit and detergents with a minimum percentage of sodium hypochlorite (or equivalent) can be used.

For a good maintenance of the equipment periodically perform the cleaning of all the external parts, by using a cloth dampened with normal neutral cleansers or simply with water.

Maintenance of the safety valve


Access the safety valve on the back side of the sterilizer.

Loosen the grained ferrule with the fingers (or with a proper tool inserted in the two holes of the ferrule) and turn counter-clockwise until you get the stop and the idle stroke.

Screw the ferrule again.

Repeat the operation at least for a couple of times.

At the end tighten the ferrule definitively.



THIS OPERATION IS NECESSARY TO GUARANTEE THE CORRECT OPERATION OF THE VALVE IN THE COURSE OF THE TIME .

AT THE END OF THE MAINTENANCE VERIFY THAT THE FERRULE IS COMPLETELY TIGHTEN AND LOCKED.

Cleaning (or replacing) the draining filter


During the use, the accumulation of residues within the filter could gradually obstruct the draining duct of the sterilization chamber.

To clean (or replace) the filter, access the rear of the sterilizer and unscrew the filter support .

Remove the filter from the support and carefully clean it under a throw of running water, using if necessary a pointed tool to remove possible extraneous bodies of greater dimensions.

If the filter cannot be recovered provide for the replacement with a new one.

Reassemble all the parts following the above procedure in opposite order.



Sterilizing the bacteriological filter

Periodically eliminate the bacterial charge of the bacteriological filter (mounted on the rear of the sterilizer) by means of a sterilization process at 121°C for porous material (program 121°C POROUS). At the end of the program screw it completely on the support.



THE STERILIZATION PROCESS DOESN'T REMOVE THE OBSTRUCTION OF THE FILTER, NEITHER PROLONGS IN ANY WAY ITS LIFE. THEREFORE PLEASE FOLLOW THE REPLACEMENT INTERVALS AS REPORTED IN THE SCHEDULED MAINTENANCE TABLE.

Replacing the bacteriological filter

At the expiration of the stated interval or whenever a visible obstruction of the filter is noticed (signaled by a color markedly tending towards the gray) unscrew the bacteriologic filter from the support, replace it with a new one that must be screwed completely on the fitting.



ONE SPARE BACTERIOLOGICAL FILTER IS SUPPLIED WITH THE EQUIPMENT. TO REQUIRE ADDITIONAL FILTERS REFER TO THE APPENDIX Z, SERVICE.

PERIODIC VALIDATION OF THE STERILIZER

As for every equipment, a decadence of the performances and components during the life, depending on the utilization type and frequency, is inevitable.

In order to guarantee a constant process safety **it is necessary to verify**, at periodic expiration (possibly yearly), the **process thermodynamic parameters** (pressure and temperature), checking if these are remaining within admitted limits or not.

The performance validation of the sterilizer is at care and **responsibility of the user**.

The reference European standards **EN 554** (*Sterilization of the medical devices - Method for the validation and systematic control of the steam sterilization*) and **EN 556** (*Sterilization of the medical devices – Requirements for the medical devices marked with “STERILE” indication*) supply an effective guide tool for carrying out the verifications on the steam sterilizers.

Since these controls require, besides a specific experience and knowledge, the use of a special tools (high accuracy sensors and probes, data logger, dedicated software, etc.), properly verified and adjusted, it is necessary to address to **companies specialized** in this activity.



“ASSISTENZA CLIENTI M.O.COM.” (SEE APPENDIX Z) IS AT YOUR DISPOSAL TO SUPPLY EVERY POSSIBLE INFORMATION CONCERNING THE PERIODIC VALIDATION OF THE STERILIZERS.

APPENDIX D - GENERAL PROBLEMS

OVERVIEW

If during the equipment use a problem or an alarm signaling occurs, it is not the case to immediately worry. This could not be caused by a breakdown, but more probably by an abnormal situation, often only transitory (for instance a black-out), or by an incorrect use.

In any case, it is important to individualize beforehand the cause of the anomaly and to carry out the opportune corrective actions autonomously or by calling the Customer Service (see Appendix Z).

To this purpose, the indications for the diagnosis and resolution of the general problems, as well as an accurate description of the alarm codes, their meaning and the consequent actions for their solution are supplied in the following.

ANALYSIS AND RESOLUTION OF THE PROBLEMS

If Your sterilizer doesn't work correctly, please perform the followings verifications before contacting the Customer Service (see Appendix Z):

PROBLEM	POSSIBLE CAUSE	PROPOSED SOLUTION
The sterilizer doesn't switch on.	The plug of the power supply cable is not correctly plugged in the socket.	Plug in correctly the power supply cable.
	Voltage lack at the output socket.	Verify the cause of the voltage lacking and try to recovery it.
	The mains switch and/or the automatic circuit breaker are in OFF state.	Set the switch in ON position.
	The mains fuses are burned.	Replace the fuses with others of the same type and value. (See Appendix A , Technical characteristics).
After the START command, the sterilization cycle doesn't start .	The equipment is performing the pre-heating phase (WARMUP).	Wait for the sterilizer reaches the correct thermal conditions for the start of the sterilization cycle. NOTE: Under normal operating conditions, the preheating average time is about 10 minutes.
The red Led MIN turns on.	The distilled water level of the internal reservoir is lower the minimum level.	Topping up the reservoir with distilled water until the lighting of the green Led MAX (or at least to turn off the red Led MIN).
The red Led Alarm turns on.	An alarm signaling has been generated with a proper code (see <i>TIMER display</i>) and message (see <i>LCD display</i>).	Verify the alarm code and operate consequently. (See following paragraphs, Alarm signaling , Alarm codes and Troubleshooting).
The safety valve intervenes.	Safety valve ferrule loosen. Anomalous over-pressure inside the chamber.	Verify the correct top ferrule locking of the safety valve. NOTE: Let the equipment cool or use gloves to avoid from burning oneself when you are touching the valve.

PROBLEM	POSSIBLE CAUSE	PROPOSED SOLUTION
At the end of the program (CYCLE COMPLETE) it is impossible to open the door.	Residual pressure inside the chamber. NOTE: LCD displays the message: NOW LEVELLING PLEASE WAIT...	Wait for some minutes until the pressure return to the 0.00 bar, and <u>retry</u> to open the door. Verify that the bacteriologic filter is not obstructed and, if necessary, replace it. The procedure for storing the environment pressure value (function SET 0 bar) has not been correctly performed. Call for a Service Intervention (see <u>Appendix Z</u>).
	The safety friction of the handle is loosen.	Push on the handle (pay attention to not damage it) against the equipment until to set it in the position parallel with the door. Retry to open the door.
	The door-locking system is remained engaged at the end of the cycle.	After having reset the alarm (code E021), verify that the handle is in the correct closing position. Call the SETUP program and the DOOR LOCK option in the menu SPECIAL for the manual command of the mechanism. (see <u>Chapter "Setting the equipment"</u>). NOTE: In case of unsuccessful, carry out the manual unlocking through the supplied tool. (see note at the end of this appendix).
Presence of water over the table under the sterilizer.	Steam leakage from the gasket.	At the end of the cycle clean with a dampened cloth the gasket and the parabola of the chamber. Check for possible damages on the gasket. Carry out a new cycle and verify the condition.
Difficulty for reaching the correct vacuum in the chamber (problems of drying, presence of water in the sterilization chamber at the end of the cycle, etc.).	Sterilization chamber draining filter obstructed.	Provide for the cleaning or replacement of the draining filter. (see <u>Appendix C "Maintenance"</u>).
	Air slots on the sides and/or on the rear obstructed or heat exchanger insufficiently aired.	Remove every possible obstruction from the air slots and from the heat exchanger. Verify the equipment is not on direct contact with walls or surfaces (see <u>Chapter "Installation"</u>).
Excessive damp on the material and/or the tools at the end of the program.	Excess of material inside the chamber.	Verify the quantity of sterilized material and make sure do not overcome the maximum admitted quantities, according to the typology of the load. (See overview table in <u>Appendix A, "Technical characteristics"</u>).
	Load not correctly arranged.	Arrange the load, particularly the wrapped one, according to the indications. (See <u>Chapter "Preparing the material to be sterilized"</u>).
	Wrong selection of the sterilization program	Choose the proper sterilization program for the type of material to be sterilized. (See Table in <u>Appendix B, "PROGRAMMES"</u>).
	Obstruction on the sterilization chamber draining filter.	Provide for cleaning or replacing the draining filter. (See <u>Appendix C, "Maintenance"</u>).

PROBLEM	POSSIBLE CAUSE	PROPOSED SOLUTION
Oxidation traces or spots on the tools	Poor quality of the tools.	Check the quality of the tools showing the problem, make sure that the material of the tools is autoclavable.
	Poor quality of the distilled water.	Empty the reservoir and fill it with high quality distilled water (possibly obtained by inverse osmosis). (See <i>Characteristics of the feeding water</i> in Appendix A , “ <i>Technical characteristics</i> ”).
	Organic or inorganic residues on the tools	Carefully clean the material before the sterilization cycle. (See Chapter “ <i>Preparing the material to be sterilized</i> ”).
	Contact among tools of different metal.	Separate the tools of different metal with proper expedients. (See Chapter “ <i>Preparing the material to be sterilized</i> ”).
	Presence of calcareous residues on the wall of the sterilization chamber and/or accessories.	Perform the cleaning of the equipment and parts as indicated. (See Appendix C , “ <i>Maintenance</i> ”).
Bluing or blackening of the materials.	Wrong selection of the sterilization program.	Verify the adequacy of the temperature of the sterilization program selected according to the material to be treated. (See <i>Table</i> in Appendix B , “ <i>PROGRAMMES</i> ”).
The printer doesn't print the report.	Wrong setup of the printer options.	Provide for setting the equipment according to the printer type used. (See Chapter “ <i>Setting the equipment</i> ”).
	No paper.	Add new paper. (see <i>printer's operating manual</i>)
	Jammed paper.	Remove the jamming. (see <i>printer's operating manual</i>)



IF THE PROBLEM PERSISTS, CALL FOR A SERVICE INTERVENTION (SEE APPENDIX Z) BY COMMUNICATING THE MODEL AND THE SERIES NUMBER OF THE STERILIZER. THESE DATA MAY BE FOUND ON THE MANUFACTURING PLATE ON THE REAR SIDE OF THE EQUIPMENT AND ON THE GUARANTEE CERTIFICATE.

APPENDIX E – ALARM SIGNALING

OVERVIEW

Each time an **anomalous condition** occurs during the sterilizer operation a warning message associated to alarm code (formed by one letter and three figures) will be generated.

The alarm codes are divided in **three categories**:

- **E = ERROR**

Wrong action and/or use, or external equipment cause.

Problem generally recoverable by the user.

Code format: **Exxx** (xxx = identification number 000 ÷ 999)

- **A = ALARM**

Failure of first level, not concerning the safety.

Problem generally recoverable on site by a specialized technician.

Code format: **Axxx** (xxx = identification number 000 ÷ 999)

- **H = HAZARD**

Failure of second level, involving the safety.

Problem generally recoverable by the Customer Service.

Code format: **Hxxx** (xxx = identification number 000 ÷ 999)

ALARM PROCEDURE



IN CASE OF ALARM, DO NOT SWITCH OFF THE STERILIZER BEFORE HAVING PERFORMED THE SYSTEM RESET (SEE PAR. "RESET OF THE SYSTEM")

The intervention of the alarm procedure causes the **interruption of the program** (or the normal operation), associated with the **alarm code** displaying, a **message** on LCD display (see table below), an **acoustic signaling** and the lighting of the **alarm red led** (intermittent).



DURING THE ALARM PROCEDURE THE DISPLAY WILL SHOW ALWAYS THE CURRENT VALUES OF THE TEMPERATURE (°C) AND PRESSURE (BAR) IN THE STERILIZATION CHAMBER

The alarm procedure is developed in order **to avoid** the user from any possible **confusion** between an anomalous cycle and a cycle correctly completed, and consequently from **unintentionally use of the not sterile material**.

The performing of the alarm procedure **differs** according to whether it occurs during or out the program execution, and is structured to drive the user up to the final **RESET** of the sterilizer.

Alarm occurring during the sterilization cycle

If the alarm occurs **during the sterilization or test program** the LCD will show:



← Depending on the alarm


← Alarm code

In case the alarm occurs in special cycle phases, the equipment enables a cleaning procedure of the hydraulic circuit, and the display will show the following message:




← Alarm code

Reached the safety conditions, the equipment enables a special procedure, and prompts the following message to manually unlock the door:



Press the key ↑ to unlock the door locking mechanism. The following message appears:



Opening the door, the user is required to reset the equipment:



Perform the **RESET** (see the next chapter), switch off the unit, analyze the problem and try for it solution.



OPENING THE DOOR AT THE END OF THE CYCLE THE REPORT (NORMAL OR EXTENDED ACCORDING TO THE TYPE OF THE OCCURRED ALARM) RELATING TO THE STERILIZATION PROGRAM INTERRUPTED AND THE ERROR SIGNALING WILL BE PRINTED. VERIFY THE PRINTOUT, SIGN IT ON THE SPECIAL SPACE AND FILE IT IN A PROPER PLACE. REFER TO THE PRINTOUT EXAMPLES IN **APPENDIX B “PROGRAMMES”**.

Alarm occurring out the sterilization cycle

If the alarm occurs out the sterilization or test program the following message will be shown:



← Alarm message

← Alarm code

Switch off the unit, analyze the problem and try for it solution.

Or, depending on the alarm type:



← Alarm message

← Alarm code

that will prompt automatically the following message:




Perform the system **RESET** (see the next chapter) , switch off the unit, analyze the problem and try for it solution.



THE ALARMS OCCURRED OUT OF THE PROGRAM WILL NOT PRODUCE ANY REPORT.

RESET OF THE SYSTEM




The **RESET** of the system can be performed in two alternative ways, according to the type of alarm occurred (see the **paragraph “List of the alarm codes”** following this appendix).

1. Hold down PROGRAM SELECTION key for about 3 seconds.
An acoustic signal will be generated as RESET confirmation;



DO NOT SWITCH OFF THE UNIT BEFORE HAVING PERFORMED THE RESET.






2. Switch off / on the unit through the mains switch.
On turning on the sterilizer the normal initial autotest will be performed.

After the RESET, and the possible technical service, the equipment will go in STAND-BY status, ready to perform a new program.

LIST OF THE ALARM CODES

The alarm codes, LCD messages and consequent RESET modes are listed in the following.

CODE	ALARM DESCRIPTION	LCD MESSAGE	RESET MODE
ERRORS (category E)			
E000	Black-out	BLACK-OUT	<div>Hold down  (> 3 seconds)</div>
E010	Door open	DOOR OPEN	
E020	Door unlocked	DOOR UNLOCKED	
E021	Door locked	DOOR LOCKED	
E030	Minimum water level in the filling tank (MIN)	WATER MIN	
E031	Maximum water level in the filling tank (MAX)	EXHAUST MAX	
E041	Too much frequent filling up of the tank (automatic filling)	FILLING PROBLEM	
E900	Unsuccessful Vacuum Test (during the LEAKAGE phase)	TEST FAILED	
E901	Unsuccessful Vacuum Test (during the WAITING phase)	TEST FAILED	
E902	Unsuccessful Vacuum Test (timeout vacuum pulses exceeded)	TEST FAILED	
E999	Manual interruption of the cycle	MANUAL STOP	
ALARMS (category A)			
A022	Problem of the door locking mechanism (OFF-OFF)	LOCKING PROBLEM	<div>Switch off the equipment</div>
A023	Problem of the door locking mechanism (ON-ON)	LOCKING PROBLEM	
A024	Problem of the door locking mechanism (ON-OFF)	LOCKING PROBLEM	
A032	Problem of the level probes	LEVEL PROBLEM	
A040	Unsuccessful reservoir filling (automatic filling)	FILLING PROBLEM	
A101	Failure of the thermo-resistor PT1 (sterilization chamber)	PTC BROKEN	
A102	Failure of the thermo-resistor PT2 (steam generator)	PTC BROKEN	
A103	Failure of the thermo-resistor PT3 (heating resistor)	PTC BROKEN	
A104	Failure of the thermo-resistor PT4 (on the wall of the sterilization chamber)	PTC BROKEN	
A111	Short-circuit of the thermo-resistor PT1 (sterilization chamber)	PTC SHORTCIRCUIT	
A112	Short-circuit of the thermo-resistor PT2 (steam generator)	PTC SHORTCIRCUIT	
A113	Short-circuit of the thermo-resistor PT3 (heating resistor)	PTC SHORTCIRCUIT	
A114	Short-circuit of the thermo-resistor PT4 (on the wall of the sterilization chamber)	PTC SHORTCIRCUIT	
A121	Instability of the thermo-resistor PT1 (sterilization chamber)	PTC FLICKERING	
A200	Pre-heating not performed within the timeout (heating resistor problem).	HEATING PROBLEM	

CODE	ALARM DESCRIPTION	LCD MESSAGE	RESET MODE
A250	1° vacuum pulse not performed within the timeout	PV1 TIMEOUT	<div>Hold down</div> <div></div> <div>(> 3 seconds))</div>
A251	1° raising towards the atmospheric pressure not performed within the timeout	ATM1 TIMEOUT	
A252	1° vacuum pulse not performed within the timeout	PP1 TIMEOUT	
A253	2° vacuum pulse not performed within the timeout	PV2 TIMEOUT	
A254	2° raising towards the atmospheric pressure not performed within the timeout	ATM2 TIMEOUT	
A255	2° vacuum pulse not performed within the timeout	PP2 TIMEOUT	
A256	3° vacuum pulse not performed within the timeout	PV3 TIMEOUT	
A257	3° raising towards the atmospheric pressure not performed within the timeout	ATM3 TIMEOUT	
A258	3° vacuum pulse not performed within the timeout	PPP TIMEOUT	
A259	PROCESS phase not started within the timeout	PROCESS TIMEOUT	
A260	Chamber leveling not performed within the timeout	PPD TIME-OUT	
HAZARDS (category H)			
H150	Failure of the pressure probe MPX	MPX BROKEN	Switch off the equipment
H160	Short-circuit/disconnection of the pressure probe MPX	MPX SHORTCIRCUIT	
H400	Ratio P_{conv}/T unbalanced ($P_{conv}>T$) (<i>PROCESS phase</i>)	P/T PROBLEM	<div>Hold down</div> <div></div> <div>(> 3 seconds)</div>
H401	Ratio T/P_{conv} unbalanced ($T>P_{conv}$) (<i>PROCESS phase</i>)	T/P PROBLEM	
H402	Temperature over the MAX limit (<i>PROCESS phase</i>)	T OVER LIMIT	
H403	Temperature under MIN limit (<i>PROCESS phase</i>)	T UNDER LIMIT	
H404	Temperature fluctuating over the limit (<i>PROCESS phase</i>)	PT1 FLUCTUATING	
H405	Pressure over the MAX limit (<i>PROCESS phase</i>)	P OVER LIMIT	
H406	Pressure under the MIN limit (<i>PROCESS phase</i>)	P UNDER LIMIT	
H 410	Process time problem (<i>PROCESS phase</i>)	TIMING PROBLEM	
H990	Over-pressure (<i>sterilization chamber, MPX</i>)	OVERPRESSURE	
H991	Over-heating (<i>sterilization chamber, PT1</i>)	OVERHEATING PT1	
H992	Over-heating (<i>steam generator, PT2</i>)	OVERHEATING PT2	
H993	Over-heating (<i>band heating resistor, PT3</i>)	OVERHEATING PT3	

ANALYSIS AND RESOLUTION OF THE PROBLEMS


Depending on the occurred alarm type, the indications to find the possible causes and restore the correct operation are listed on the following table:

CODE	POSSIBLE CAUSE	SUGGESTED SOLUTION
ERRORS (category E)		
E000	Sudden interruption of the electric power supply (black-out).	Wait for the return of the electric power supply and perform the RESET according to the instructions. Verify the correct sterilization of the load before using the treated material.
	Accidental switching off of the main switch and/or disconnection of the plug from the a.c. socket.	Reconnect the plug and/or switch on the equipment, and perform the RESET according to the instructions. Verify the correct sterilization of the load before using the treated material.
	Mains fuses burned.	Replace the fuses of same type and value (see <i>Table in Appendix A, "Technical characteristics"</i>). Switch on the equipment and perform the RESET according to the instructions. Verify the correct sterilization of the load before using the treated material.
E010	Door open (or not correctly close) on the start of the cycle (START).	Perform the RESET according to the instructions. <u>Correctly</u> close the door and restart the cycle.
	Failure of the door micro-switch.	Call for a Service Intervention (see Appendix Z).
E020	Door locking mechanism not correctly activated on the start of the cycle (START).	Perform the RESET according to the instructions. <u>Correctly</u> close the door, pushing completely the handle, and try to restart again the program
	Failure of the door locking mechanism.	Call for a Service Intervention (see Appendix Z).
E021	Door locking system not released at the end of the cycle (CYCLE COMPLETE).	Perform the RESET according to the instructions. Call the SETUP program and select the DOOR LOCK option in the menu SPECIAL in order to operate manually the mechanism. (see <i>Chapter "Setting the equipment"</i>). NOTE: In case of unsuccessful result, manually perform the releasing of the mechanism through the supplied tool.
	Failure of the door locking mechanism.	Call for a Service Intervention (see Appendix Z).
E030	Internal reservoir water level under the minimum (MIN).	Perform the RESET according to the instructions. Provide for filling up the water until the turning on of the Led MAX (or turning off at least of the Led MIN).
	Failure of the MIN level signaling.	Call for a Service Intervention (see Appendix Z).
E031	Draining tank water level over the MAX level.	Perform the RESET according to the instructions. Provide for draining the tank, letting a minimum level of water as marked to remain.
	Failure of the MAX level signaling.	Call for a Service Intervention (see Appendix Z).

CODE	POSSIBLE CAUSE	SUGGESTED SOLUTION
E041	Connecting pipe between external tank and sterilizer not correctly connected.	Perform the RESET according to the instructions. Verify the connection of the water filling pipe, tighten both ends on the fitting through the supplied plastic clips. Eliminate possible obstructions on the pipe route.
	Failure of the water pump.	Call for a Service Intervention (see Appendix Z).
	Problem on the hydraulic circuit.	Call for a Service Intervention (see Appendix Z).
E900	Air leakage on the gasket	Perform the RESET according to the instructions. Carefully clean the gasket with a clean cloth of cotton dampened with water. Restart the program.
	Sterilization chamber too much warm.	Perform the RESET according to the instructions. Carry out the test with the sterilization chamber at moderate temperature (for instance in the morning, at the first switching on of the equipment).
	Problem on the hydraulic circuit.	Call for a Service Intervention (see Appendix Z).
E901	Excess of damp inside the sterilization chamber.	Perform the RESET according to the instructions. Carefully dry the sterilization chamber and restart the program.
	Air leakage on the gasket	Perform the RESET according to the instructions. Carefully clean the gasket with a clean cloth of cotton dampened with water. Restart the program.
	Problem on the hydraulic circuit.	Call for a Service Intervention (see Appendix Z).
E902	Excess of damp in the sterilization chamber.	Perform the RESET according to the instructions. Carefully dry the inside of the sterilization chamber and again restart the program.
	Air leakage on the gasket	Perform the RESET according to the instructions. Carefully clean the gasket with a clean cloth of cotton dampened with water. Restart the program.
	Failure of the vacuum pump.	Call for a Service Intervention (see Appendix Z).
	Problem on the hydraulic circuit.	
E999	Manual interruption of the sterilization or test program. (see also Chapter “Running the sterilization program”)	Perform the RESET according to the instructions. Verify the correct load sterilization before using the treated material.
ALARMS (category A)		
A022	Failure of the door locking mechanism	Call for a Service Intervention (see Appendix Z).
A023	Failure of the door locking mechanism	
A024	Failure of the door locking mechanism	

CODE	POSSIBLE CAUSE	SUGGESTED SOLUTION
A032	Water level signaling connector not plugged in. Failure of the level signaling(s)	Call for a Service Intervention (see Appendix Z).
A040	Lacking of water in the external tank (automatic filling).	Perform the RESET according to the instructions. Fill the tank with sufficient distilled water, and remember to periodically verify the level.
	Connecting pipe between the external tank and sterilizer not correctly connected.	Perform the RESET according to the instructions. Verify the connection of the water filling pipe, tighten both ends on the fitting through the supplied plastic clips. Eliminate possible obstructions on the pipe route.
	Failure of water pump.	
A101	Failure of the temperature probe of the sterilization chamber (PT1).	Call for a Service Intervention (see Appendix Z).
A102	Failure of the temperature probe of the steam generator (PT2).	
A103	Failure of the temperature probe of the heating resistor (PT3).	
A104	Failure of the temperature probe of the chamber wall (PT4).	
A111	Wrong connection of the temperature probe (sterilization chamber) to the pcb connector. Short-circuit of the temperature probe (sterilization chamber).	
A112	Wrong connection of the temperature probe (steam generator) to the pcb connector. Short-circuit of the temperature probe (steam generator).	
A113	Wrong connection of the temperature probe (heating resistor) to the pcb connector. Short-circuit of the temperature probe (heating resistor).	
A114	Wrong connection of the temperature probe (chamber wall) to the pcb connector. Short-circuit of the temperature probe (chamber wall).	
A121	Unstable operation of the temperature probe (sterilization chamber).	

CODE	POSSIBLE CAUSE	SUGGESTED SOLUTION
A200	Release of the steam generator safety thermostat.	<p>Provide for the manual reset of the thermostat installed under the equipment (see Chapter “Product description”).</p> <p>Unscrew the protection plastic cap, push on the red button for a click and replace the cap.</p> <p>Switch off (RESET) and switch on again the equipment.</p>
	Release of the heating resistor safety thermostat. Wrong operation of the steam generator or heating resistor.	<p>Call for a Service Intervention (see Appendix Z).</p>
A250	Presence of water or condensation inside the sterilization chamber.	<p>Perform the RESET according to the instructions.</p> <p>Carefully dry the sterilization chamber and restart the program.</p> <p>Do not introduce material soaked by water or liquid into the chamber.</p>
	Obstruction on the sterilization chamber draining filter.	<p>Provide for cleaning or replacing the draining filter.</p> <p>(See Appendix C, “Maintenance”).</p>
	Air leakage on the gasket	<p>Perform the RESET according to the instructions.</p> <p>Carefully clean the gasket with a clean cloth of cotton dampened with water.</p> <p>Restart the program.</p>
	Failure of the vacuum pump.	<p>Call for a Service Intervention (see Appendix Z).</p>
	Problem on the hydraulic circuit.	
A251	Wrong operation of the water pump. Problem on the hydraulic circuit.	
	Release of the steam generator safety thermostat.	<p>Provide for the manual reset of the thermostat installed under the equipment (see Chapter “Product description”).</p> <p>Unscrew the protection plastic cap, push on the red button until a click and replace the cap.</p> <p>Switch off (RESET) and switch on the equipment.</p>
	Release of the heating resistor safety thermostat. Wrong operation of the steam generator or heating resistor.	<p>Call for a Service Intervention (see Appendix Z).</p>
A252	Air leakage on the gasket	<p>Perform the RESET according to the instructions.</p> <p>Carefully clean the gasket with a clean cloth of cotton dampened with water.</p> <p>Restart the program.</p>
	Excess of load.	<p>Perform the RESET according to the instructions.</p> <p>Verify the quantity of material into the sterilization chamber and pay attention to not exceed the maximum admitted quantity according to the load typology.</p> <p>(see Table in Appendix A, Technical characteristics).</p>
	Problem of the hydraulic circuit.	<p>Call for a Service Intervention (see Appendix Z).</p>




CODE	POSSIBLE CAUSE	SUGGESTED SOLUTION
A252 (continue)	Release of the steam generator safety thermostat.	Provide for the manual reset of the thermostat installed under the equipment (see Chapter “Product description”). Unscrew the protection plastic cap, push on the red button until a click and replace the cap. Switch off (RESET) and switch on again the equipment.
	Release of the heating resistor safety thermostat. Wrong operation of the steam generator or heating resistor.	Call for a Service Intervention (see <u>Appendix Z</u>).
A253	Presence of water or condensation into the sterilization chamber.	Perform the RESET according to the instructions. Carefully dry the sterilization chamber and restart the program. Do not introduce material soaked with water or liquid into the chamber.
	Air leakage on the gasket	Perform the RESET according to the instructions. Carefully clean the gasket with a clean cloth of cotton dampened with water. Restart the program.
	Failure of the vacuum pump.	Call for a Service Intervention (see <u>Appendix Z</u>).
	Problem on the hydraulic circuit.	
A254	Wrong operation of the water pump. Problem on the hydraulic circuit.	Call for a Service Intervention (see <u>Appendix Z</u>).
	Release of the steam generator safety thermostat.	Provide for the manual reset of the thermostat installed under the equipment (see Chapter “Product description”). Unscrew the protection plastic cap, push on the red button until a click and replace the cap. Switch off (RESET) and switch on again the equipment.
	Release of the heating resistor safety thermostat. Wrong operation of the steam generator or heating resistor.	Call for a Service Intervention (see <u>Appendix Z</u>).
A255	Air leakage on the gasket	Perform the RESET according to the instructions. Carefully clean the gasket with a clean cloth of cotton dampened with water. Restart the program.
	Excess of load.	Perform the RESET according to the instructions. Verify the quantity of material in the sterilization chamber and pay attention to not exceed the maximum admitted quantity according to the load typology. (see <i>Table in <u>Appendix A, Technical characteristics</u></i>).
	Problem of the hydraulic circuit.	Call for a Service Intervention (see <u>Appendix Z</u>).
	Release of the steam generator safety thermostat.	Provide for the manual reset of the thermostat installed under the equipment (see Chapter “Product description”). Unscrew the protection plastic cap, push on the red button until a click and replace the cap. Switch off (RESET) and switch on again the equipment.

CODE	POSSIBLE CAUSE	SUGGESTED SOLUTION
A255 (continue)	Release of the heating resistor safety thermostat.	Call for a Service Intervention (see Appendix Z).
	Wrong operation of the steam generator or heating resistor.	
A256	Presence of water or condensation in the sterilization chamber.	Perform the RESET according to the instructions. Carefully dry the sterilization chamber and restart the program. Do not introduce material soaked with water or liquid into the chamber.
	Air leakage on the gasket	Perform the RESET according to the instructions. Carefully clean the gasket with a clean cloth of cotton dampened with water. Again start the program.
	Failure of the vacuum pump.	Call for a Service Intervention (see Appendix Z).
	Problem on the hydraulic circuit.	
A257	Wrong operation of the water pump. Problem on the hydraulic circuit.	Call for a Service Intervention (see Appendix Z).
	Release of the steam generator safety thermostat.	Provide for the manual reset of the thermostat installed under the equipment (see Chapter “Product description”). Unscrew the protection plastic cap, push on the red button until a click and replace the cap. Switch off (RESET) and switch on again the equipment.
	Release of the heating resistor safety thermostat. Wrong operation of the steam generator or heating resistor.	Call for a Service Intervention (see Appendix Z).
A258	Air leakage on the gasket	Perform the RESET according to the instructions. Carefully clean the gasket with a clean cloth of cotton dampened with water. Restart the program.
	Excess of load.	Perform the RESET according to the instructions. Verify the quantity of material in the sterilization chamber and pay attention to not exceed the maximum admitted quantity according to the load typology. (see <i>Table in Appendix A, Technical characteristics</i>).
	Problem on the hydraulic circuit.	Call for a Service Intervention (see Appendix Z).
	Release of the steam generator safety thermostat.	Provide for the manual reset of the thermostat installed under the equipment (see Chapter “Product description”). Unscrew the protection plastic cap, push on the red button until a click and replace the cap. Switch off (RESET) and switch on again the equipment.
	Release of the heating resistor safety thermostat. Wrong operation of the steam generator or heating resistor.	Call for a Service Intervention (see Appendix Z).


CODE	POSSIBLE CAUSE	SUGGESTED SOLUTION
A259	Excess of load.	Perform the RESET according to the instructions. Verify the quantity of material in the sterilization chamber and pay attention to not exceed the maximum admitted quantity according to the load typology. (see <i>Table in Appendix A, Technical characteristics</i>).
	Air leakage from the gasket	Perform the RESET according to the instructions. Carefully clean the gasket with a clean cloth of cotton dampened with water. Restart the program.
	Problem on the hydraulic circuit.	Call for a Service Intervention (see <u>Appendix Z</u>).
A260	Problem on the hydraulic circuit	
HAZARDS (category H)		
H150	Break of the pressure probe (MPX).	Call for a Service Intervention (see <u>Appendix Z</u>).
H160	Wrong connection of the pressure probe cable (MPX) to the p.c.b. connector. Short-circuit of the pressure probe (MPX).	
H400	Problem of the hydraulic circuit.	
H401	Problem of the hydraulic circuit.	
H402	Wrong operation of the steam generator. Problem of the hydraulic circuit.	
H403	Wrong operation of the steam generator. Problem of the hydraulic circuit.	
H404	Wrong operation of the steam generator. Problem of the hydraulic circuit.	
H405	Wrong operation of the steam generator. Problem of the hydraulic circuit.	
H406	Wrong operation of the steam generator. Problem of the hydraulic circuit.	
H410	Process timing problem	
H990	General operation problem.	
H991	General operation problem.	
H992	General operation problem.	
H993	General operation problem.	

APPENDIX F – DRAWINGS

ELECTRICAL DRAWING



HYDRAULIC DRAWING



APPENDIX G – DECLARATION OF CONFORMITY

Sistemi innovativi di sterilizzazione

DECLARATION OF CONFORMITY

Application of the EEC Directives 93/42 - 89/336 - 73/23

Name of the Manufacturer: **M.O.COM. S.r.l. - Manifattura Odontoiatrica Complementare**Address of the Manufacturer: **Via delle Azalee, 1 - 20090 Buccinasco (MI) - ITALY**Product description: **Steam Sterilizer**Model: **millennium[®] B_μ**Made in: **ITALY**

The undersigned hereby declares that the above mentioned goods entirely

conform

to the EEC Directives 93/42 - 89/336 - 73/23 (and following updates).

Reference standards: EN 61010-1 EN 61010-1-A2 EN 61010-2-041
CEI EN 50081-1 CEI EN 50082-1
EN 55014 EN 55022 EN 60555-2 (CEI 77-3) EN 60555-3 (CEI 77-4)
EN 61000-4-2 EN 61000-4-3 EN 61000-4-4
ASME VIII Div. 1 (Add. 1999) DIN 58946 T5 TRD 421 TRD 511
prEN 13060: 2004

31/03/2004

Date

Signature

Alfio VILLA
Name and SurnameLegal Agent
Function

APPENDIX H –OPERATOR’S NOTES

[illegible]

APPENDIX Z – CUSTOMER SERVICE

**FOR ANY SERVICE INTERVENTION
BOTH DURING AND OUT THE PRODUCT
WARRANTY PERIOD
PLEASE CALL DIRECTLY**

THE CUSTOMER SERVICE

**OF THE AGENT OR RETAILER
THAT SUPPLIED THE EQUIPMENT**

M.O.COM. Ltd. Co. is at complete disposal of the Customers for any technical additional information concerning the product, as well as for any suggestions on the steam sterilization procedures.

Please contact at the following address:

M.O.COM. Srl
Assistenza Clienti
Via delle Azalee, 1
20090 Buccinasco (MI)
ITALY

Tel. (+39) 02-45701505
Fax (+39) 02-45701258
e-mail at@mocom.it

To help us in product quality and service improving, please send Your comments and/or suggestions to the followings **e-mail** addresses:

uc@mocom.it (marketing and sales)

Besides, You can send comments and/or suggestions by **letter or fax** to the above indicated address.

Thanks in advance for the valuable assistance that you would supply to us.